MANUEL D'UTILISATION • USER GUIDE • GEBRAUCHSANLEITUNG • MANUAL DE USO • ISTRUZIONI D'USO • INSTRUÇÕES DE USO • BRUKSANVISNING
GEBRUIKSAANWIJZING • INSTRUKCJA OBSŁUGI • ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ • ИНСТРУКЦИЯ ДЛЯ ПОЛЬЗОВАТЕЛЯ • UŽIVATELSKÝ MANUÁL • 使用说明 • 使用説明書

Respireo **SOFT**



Nasal

Masque nasal avec fuites calibrées • Nasal mask vented • Nasenmaske mit kalibrierten luftauslassöffnungen • Mascarilla nasal con
orificios de exhalación calibrados • Maschera nasale con fori di esalazione calibrati • Máscara nasal com furos de exalação calibrados
Näsmask med kalibrerade utandningshål • Neusmasker met gekalibreerde ventilatieopeningen • Maska nosowa z kalibrowanymi
otworami wentylacyjnymi • Ρινική μάσκα με ρυθμισμένες οπές εξαερισμού • Носовая маска с калиброванными выдыхательными
отверстиями • Nosní maska s kalibrovanými vydechovacími otvory • 带有标准通气口的 鼻面罩 • 鼻マスク 通気穴付き


FRANÇAIS.....	p. 8
ENGLISH.....	p. 12
DEUTSCH.....	p. 15
ESPAÑOL.....	p. 20
ITALIANO.....	p. 24
PORTUGUÊS.....	p. 29
SVENSKA.....	p. 33
NEDERLANDS.....	p. 37
POLSKI.....	p. 41
ΕΛΛΗΝΙΚΑ.....	p. 45
РУССКИЙ.....	p. 50
ČEŠTINA.....	p. 54
中文.....	p. 57
日本語.....	p. 61

A


B


1

2


3


B


④











⑤

⑥

C**1****2****D****1****3****2****4**

E**1****2****3****4**

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1 – INTENDED USE

Respireo SOFT Nasal is a nasal mask with calibrated vent ports (Vented) intended for single-patient reuse in the home environment and multipatient reuse in the hospital environment .

Respireo SOFT Nasal is to be used by adult patients (weighing >30 kg) for whom a continuous positive pressure therapy (CPP) or a respiratory support with positive expiratory pressure (PEP) have been prescribed.

⚠ 2 - WARNINGS

- 2.1 The mask must only be used with positive pressure generators recommended by a physician or respiratory therapist.
- 2.2 The mask must only be used with the appliance turned on. The flow through the vent ports is guaranteed if the appliance provide at least a pressure of 4 cmH₂O. With low pressures, the flow through the vent ports may be insufficient, which could lead to partial re-inhalation.
- 2.3 The vents must be kept clear.
- 2.4 If you experience discomfort or irritation, or suffer from an allergic reaction to any of the components of the mask, please report this

to your doctor or respiratory therapist.

- 2.5 If oxygen is used, do not smoke or use open flame.
- 2.6 If oxygen is used, the flow must be turned off when the device is not operating, in order to prevent an accumulation of unused oxygen and a risk of fire.
- 2.7 Do not use the mask in the presence of nausea or vomiting.
- 2.8 The mask must be kept away from light.
- 2.9 The mask must be cleaned and/or disinfected. For further details refer to the section "Cleaning and disinfecting the mask".
- 2.10 Check the mask before use. If there is any visible damage, report it to the dealer.
- 2.11 Replace the mask if it shows any signs of deterioration (e.g. cracks, scratches, fractures or tears).
- 2.12 The mask must only be used in accordance with these instructions. No other instructions are applicable.
- 2.13 Some of the components of the mask are small enough to be swallowed by children and must never be left unattended. An adult must always be present.

3 - COMPONENTS OF THE MASK (Fig. A)

- a) cushion; b) mask frame; c) elbow; d) swivel connector; e) headgear; f) buckle.

4 – FITTING (Fig.B)

- 4.1 For the first installation, first of all remove the packaging then attach the cushion (a) to the mask frame (b) keeping the cushion triangle aligned with the mask frame triangle. To properly mount the cushion (a), put the fingers and push as shown in the picture (Fig.B1).
Note: Be sure that the lines on the mask frame (b) and on the cushion (a) are aligned; if not, the cushion is not mounted properly. Please mount the cushion again.
- 4.2 Keep the mask in one hand and align the bottom part of the cushion opening to the bottom part of the nose then with the other

- hand pull the headgear (e) over and behind the head (Fig.B2).
- 4.3 Connect the two lower straps of the headgear (e) to the hooks. Ensure the lower straps are under ears (Fig.B3).
 - 4.4 Gently tighten the straps to ensure stability of the mask, first the strap on the head, second the upper straps (on the forehead) and third the lower straps (Fig.B4).

Notes:

- Repeat steps if necessary
 - The lower straps are intended to hold the mask in position. To ensure the stability of the mask, they must not be kept too loose.
- 4.5 Connect the swivel connector (d) to the device tubing and then switch on the device at the prescribed pressure (Fig.B5).
 - 4.6 Lye down in resting position and readjust the forehead straps (Fig.B6).

5 – REMOVING THE MASK (Fig.C)

- 5.1 Disconnect both lower straps of the headgear (e) from the mask frame hooks (Fig.C1).
Note: Keep the hook tabs in position in order to keep your settings for the next time.
- 5.2 Pull the mask up from the nose and lift the headgear (e) over the head (Fig.C2).

6 – DISASSEMBLING THE MASK (Fig.D)

Note: Before disassembling the mask, disconnect the patient circuit from the mask.

- 6.1 Disconnect the straps of the headgear (e) from the mask frame (b) (Fig.D1).
Note: Do not undo the hook tabs in order to keep the setting for the next time.
- 6.2 Separate the cushion (a) from the mask frame (b) (Fig.D2).
- 6.3 Disconnect the elbow (c) from the mask frame (b) (Fig.D3).
- 6.4 Disconnect the swivel connector (d) from the elbow (c) (Fig.D4).

7 - CLEANING THE MASK

7.1 – CLEANING AT HOME FOR REUSE BY THE SAME PATIENT
First, disassemble all the components (see DISASSEMBLING THE MASK).

After each use:

Carefully clean the cushion (a) and the mask frame (b) in warm water (about 30°C) using mild soap. Thoroughly rinse under the tap and leave in the air to dry, out of direct sunlight.

Every 7 days:

Carefully clean the cushion (a), the mask frame (b) and the headgear (e) in warm water (about 30°C) using mild soap. Thoroughly rinse under the tap and leave in the air to dry, out of direct sunlight.

Notes:

- Do not clean the mask with solutions containing alcohol, aromatic compounds, humidifiers, antibacterial agents, hydrating agents or aromatic oils as they would damage the product and reduce its life.
- Do not leave the mask immersed for more than 10 minutes.
- Do not clean the mask in the dishwasher.
- Check carefully all the components to make sure there are no soap deposits or foreign bodies inside the mask.

7.2 – REPROCESSING THE MASK FOR RE-USE WITH OTHER PATIENTS

Notes:

- The mask must always be reprocessed before being used by another patient.
- The headgear (e) must be replaced each time.

The mask can be regenerated up to 20 times, using one of the following methods:

- high-level thermal disinfection;
- high-level chemical disinfection;

Full details of the mask regeneration process can be obtained by

logging on to www.device.airliquidehealthcare.com or by contacting your dealer.

8 - REASSEMBLY THE MASK (Fig.E)

- 8.1 Insert the swivel connector (d) in the elbow (c) (Fig.E1).
- 8.2 Insert the elbow (c) in the mask frame (b) (Fig.E2).
- 8.3 If the buckle (f) has been detached from the headgear, insert the straps located over the head into it by keeping the light color side of the straps facing out.
- 8.4 Insert the two upper straps into the slots located on the forehead support of the mask frame (b) (Fig.E3).
Note: Do not open the headgear hooks in order to do not readjust the straps.
- 8.5 Attach the cushion (a) to the mask frame (b) keeping the cushion triangle aligned with the mask frame triangle (Fig.E4).
Note: Be sure that the lines on the mask frame (b) and on the cushion (a) are aligned; if not, the cushion is not mounted properly. Please mount the cushion again.

9 - TROUBLESHOOTING

PROBLEM	POSSIBLE CAUSES	SOLUTIONS
• No air comes out of the vent ports	• Flow generator is off or not working properly	• Switch on the flow generator • Regulate the flow generator
• Difficulty of breathing	• Patient circuit tube is not well inserted	• Connect the tube to the flow generator and mask.
• The mask is noisy	• Vent ports are or partially obstructed	• Clean the vent ports with a soft bristle brush
	• Mask is not properly assembled	• Disassemble the mask, then reassemble according to the instruction
	• Mask is not correctly positioned	• Refit or reposition the mask and set the headgear according the instruction of fitting

PROBLEM	POSSIBLE CAUSES	SOLUTIONS
• Mask leaks around the face	• Mask is not positioned correctly • Headgear straps are not correctly set (over tighten or too loose).	• Refit or reposition the mask and set the headgear according the instruction of fitting


10 - DISPOSAL INSTRUCTIONS

The components of the mask are made of the following materials and must be disposed of in accordance with the applicable laws.

COMPONENT	MATERIAL
Cushion (a)	Silicone + Polybutylene
Mask frame (b)	Polycarbonate + Silicone
Elbow (c)	Polycarbonate
Swivel connector(d)	Nylon
Headgear (e)	Nylon+Poliurethane
Buckle (f)	Nylon

11 - TECHNICAL SPECIFICATIONS

PRESSURE-FLOW CURVE

Note: Due to manufacturing differences, the exhalation flow rate may change. (±10%).

Normative	CE marking in accordance with European Directive 93/42/EEC, class IIa and subsequent amendments EN ISO 14971, ISO 5356-1, ISO 10993-1, ISO 17510, EN ISO 17664
Operating pressure	4 cmH ₂ O – 40 cmH ₂ O
Resistance	At 50 lpm: 0,24 cmH ₂ O – at 100 lpm 0,93 cmH ₂ O
Dead space	98 ml (Large cushion)
A-weighted sound pressure level according to ISO 3744	16,7 dB(A) ± 3,2 dB(A)
A-weighted sound power level according to ISO 3744	24,7 dB(A) ± 3,2 dB(A)
Operating temperature	5°C – 40°C
Operating humidity	15 % – 95 %
Gross dimensions	135 mm x 120 mm x 80 mm

12 - WARRANTY

This product by Air Liquide Medical Systems S.r.l. is guaranteed against material and manufacturing defects for a period of 120 days as from the date of purchase, provided it is used in accordance with the operating instructions and for the intended purpose only. If the product fails under conditions of normal use, Air Liquide Medical Systems S.r.l. will repair or replace, at its option the defective product or any of its components.

The manufacturer is liable for the safety, reliability and performance of the apparatus, provided that it is used in accordance with the instructions and for the intended purpose only, and that any repairs are made by the manufacturer or an approved service centre. Throughout the warranty period keep the original packing. Warranty claims on effective product must be made by the consumer

contacting the dealer.
For further information on your warranty rights contact your dealer.

COPYRIGHT

The information contained in this booklet may not be used for anything other than the intended purpose. This booklet is the property of Air Liquide Medical Systems S.r.l. and no part of it may be reproduced without the manufacturer's written permission. All rights reserved.

TECHNICAL UPDATES

Air Liquide Medical Systems S.r.l. periodically reviews and modifies all its medical devices in order to improve their performance, safety and reliability. The instruction booklets are updated to ensure that they continue to be consistent with the characteristics of the devices marketed. If the booklet accompanying the apparatus is mislaid, a replacement copy can be obtained from the Manufacturer by quoting the data shown on the label.

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AIR LIQUIDE MEDICAL SYSTEMS S.r.l.