MANUEL D'UTILISATION • USER GUIDE • GEBRAUCHSANLEITUNG • MANUAL DE USO • ISTRUZIONI D'USO • INSTRUÇÕES DE USO • BRUKSANVISNING
GEBRUIKSAANWIJZING • INSTRUKCJA OBSŁUGI • ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ • ИНСТРУКЦИЯ ДЛЯ ПОЛЬЗОВАТЕЛЯ • UŽIVATELSKÝ MANUÁL • 使用说明 • 使用説明書


Respireo Primo F

MASQUE BUCCO-NASAL • FACIAL MASK • MUND-NASEN-MASKE • MASCARILLA ORO-NASAL • MASCHERA ORO-NASALE
MÁSCARA ORO-NASAL • HELMASK • MOND-NEUSMASKER • MASKA USTNO-NOSOWA • ΠΙΝΙΚΗ ΣΤΟΜΑΤΙΚΗ ΜΑΣΚΑ
ПОТО-НОЦОВАЯ МАСКА • NOSO-ÚSTNÍ MASKA • 口鼻面罩 • 鼻口腔マスク


FRANÇAIS.....	p. 7
ENGLISH.....	p. 11
DEUTSCH.....	p. 15
ESPAÑOL.....	p. 20
ITALIANO.....	p. 24
PORTUGUÊS.....	p. 29
SVENSKA.....	p. 33
NEDERLANDS.....	p. 37
POLSKI.....	p. 42
ΕΛΛΗΝΙΚΑ.....	p. 46
РУССКИЙ.....	p. 51
ČEŠTINA.....	p. 56
中文.....	p. 60
日本語.....	p. 62

A


B


1


2





3


4


B

C

D

1 - UTILISATION

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1 - APPLICATION

Respireo Primo F is a mouth-and-nose mask with calibrated vent ports for use by a single patient.

Respireo Primo F is designed for use by patients who have been prescribed continuous positive airway pressure (CPAP) or Bi-level positive airway pressure (BiPAP) treatment at home or a hospital.

Respireo Primo F is suitable for adult patients (>30Kg).

⚠ 2 - IMPORTANT NOTES

- 2.1 The mask must only be used with CPAP or Bi-level flow generators recommended by a doctor or respiratory therapist.
- 2.2 The operating pressure range is 4-40 cmH₂O.
- 2.3 The mask must only be used with the appliance on and operating at a pressure of at least 4 cmH₂O. With low pressures using a CPAP appliance, the flow through the vent ports may be insufficient to eliminate all the gas exhaled from the tube, which could lead to partial re-inhalation.
- 2.4 The vents in the mask must never be obstructed.
- 2.5 If you experience discomfort or irritation, or suffer an allergic reaction to any of the components of the mask, please report this to your doctor or respiratory therapist.
- 2.6 If oxygen is used, the supply must be turned off when the appliance is not in operation, to prevent oxygen accumulating inside, which is a fire risk.
- 2.7 Do not smoke or light naked flames when using oxygen.
- 2.8 Do not use the mask in the presence of nausea or vomiting.
- 2.9 The mask must be cleaned and/or disinfected. For further details refer to the section 7 - CLEANING AND DISINFECTING THE MASK.
- 2.10 Check the mask before use. If there is any visible damage, report it to the dealer.
- 2.11 The mask must only be used in accordance with these instructions. No other instructions are applicable.

2.12 Some of the components of the mask are small enough to be swallowed by children and must never be left unattended. An adult must always be present.

3 - COMPONENTS OF THE MASK (Fig. A)

a) mouth-and-nose cushion; b) main body; c) forehead support; d) forehead rest e) adjusting knob; f) adjusting screw; g) cap; h) anti-suffocation fitting; i) swivel fitting; l) headgear; m) clip.

4 - HOW TO PUT ON THE MASK (Fig. B)

The mask comes fully assembled ready for use. Proceed as follows to put on the mask.

Turn the adjusting knob (e) until the forehead support (c) is fully extended (Fig. 1).

Hold the mask firmly in one hand and extend the headgear (l) with the other enough to pull it over your head (Fig. 2).

Place the mask gently on your face so that the cushion (a) fits snugly round the mouth and nose, and position the headgear (l) on your head (Fig. 3).

Clip the two clips (m) onto the main body (b) (Fig. 4).

Tighten the upper straps first, then the lower straps, until the mask fits on your face comfortably, without over tightening (Fig. 5).

Connect the mask to a CPAP or a Bi-level appliance at the prescribed pressure (Fig. 6)

After lying down, turn the adjusting knob (e) until the mask feels comfortable but still provides a good seal (Fig. 7).

Readjust the lower straps if necessary to eliminate any leaks, but make sure the mask is not too tight.

5 - HOW TO REMOVE THE MASK (Fig. C)

To remove the mask, press the button on the front (Fig. C) to release one of the clips (m). Lift the headgear over your head and remove the mask.

The other straps do not need to be loosened or disconnected.

6 - HOW TO DISASSEMBLE THE MASK

Before disassembling the mask, disconnect the patient circuit from the CPAP or Bi-level appliance.

Make a note of the position of the adjusting knob (e) for reference next time.

Unclip the two clips (m) from the main body (b).

Undo the four straps and detach the headgear (l) from the main body (b) and the clips (m).

Pull the anti-suffocation fitting (h) off the main body (b).

Pull the swivel fitting (i) off the anti-suffocation fitting (h).

Remove the cushion (a) from the main body (b). Remove the cap (g) from the main body (b).

Detach the forehead rest (d) from the forehead support (c).

7 - CLEANING AND DISINFECTING THE MASK

The mask must be cleaned and/or disinfected each time before use.

N.B. The headgear should be cleaned once a week.

CLEANING THE MASK

First wash your hands thoroughly.

Disassemble all the components (see 6 - HOW TO DISASSEMBLE THE MASK) and wash carefully in water at about 30°C using mild soap.

Rinse thoroughly with water and leave in the air to dry.

The mask can be cleaned in a washing machine at 30°C (programme for delicate fabrics) but must not be spin- or tumble-dried. The mask can be washed in this way up to 20 times.

DISINFECTING THE MASK

The mask must be taken apart for disinfection (see 6 - HOW TO DISASSEMBLE THE MASK). Immerse the components in a cold solution of disinfectant (e.g. a solution containing up to 2% sodium hypochlorite), following the manufacturer's instructions.

IMPORTANT. Check that the disinfectant solution is compatible with the materials of which the mask is made (see 12 - DISPOSAL

INSTRUCTIONS).

N.B. The headgear (l) must not be disinfected.

8 - PRECAUTIONS

Do not clean the mask with solutions containing alcohol, aromatic compounds, humidifiers, antibacterial agents or aromatic oils. Such solutions would damage the product and reduce its life. Prolonged exposure to direct sunlight or heat may cause damage to the mask.

If any visible damage is noted (cracks, tears, etc.), the damaged part must be replaced.

Check carefully to make sure there are no soap deposits or foreign bodies inside the mask.

Check regularly that the vents are not obstructed. Clean with a soft brush if necessary.

9 - HOW TO REASSEMBLE THE MASK

Insert the anti-suffocation fitting (h) in the main body (b) of the mask, pressing gently.

Remount the forehead rest (d) on the forehead support (c), keeping the curved edge facing upwards.

Insert the swivel fitting (i) in the anti-suffocation fitting (h), pressing gently.

Fit the grooved side of the cushion (a) into the frame of the main body (b). Fit on the top first, then the sides of the mask, pressing the cushion in until the whole edge is up against the six stops on the main body (Fig. D).

Fit on the cap (g).

Pass the two lower straps of the headgear (l) through the slots in the clips (m).

Attach the two upper straps of the headgear (l) to the main body (b). Clip the two clips (m) back onto the main body (b).

Turn the adjusting knob (e) to the reference position you noted previously.

10 - PRESSURE/OXYGEN PORT

If additional oxygen or a pressure reading is required, use the port in the main body (b).

Open the cap (g) by lifting upwards, without removing it.

Push the pressure/oxygen tube firmly into the port.

11 - TROUBLE-SHOOTING CHART

PROBLEM	POSSIBLE CAUSES	SOLUTIONS
• No air or insufficient air comes out of the mask	<ul style="list-style-type: none"> • Flow generator off • Flow generator regulated wrongly • Patient circuit tube not inserted 	<ul style="list-style-type: none"> • Switch on the flow generator • Regulate the flow generator • Connect the tube to the flow generator and mask
• The mask is noisy	<ul style="list-style-type: none"> • Cushion mounted incorrectly and air leak • Anti-suffocation fitting dirty or damaged 	<ul style="list-style-type: none"> • Mount the cushion correctly • Clean or replace the anti-suffocation fitting
• The mask does not fit the face snugly	<ul style="list-style-type: none"> • Mask applied incorrectly 	<ul style="list-style-type: none"> • Reposition the headgear
• Exhaled air does not flow out and difficulty inhaling	<ul style="list-style-type: none"> • Vent ports obstructed 	<ul style="list-style-type: none"> • Remove the obstruction


12 - DISPOSAL INSTRUCTIONS

The components of the mask are made of the following materials and must be disposed of in accordance with the applicable laws.

COMPONENT	MATERIAL
Main body	Polycarbonate
Forehead support	Polycarbonate
Adjusting knob	Polypropylene
Adjusting screw	Polycarbonate
Anti-suffocation fitting	Polycarbonate
Anti-suffocation valve	Silicone
Swivel fitting	Polycarbonate
Forehead rest	Silicone
Cushion	Silicone
Cap	Silicone
Headgear	Nylon/polyurethane
Clip	Nylon

13 - TECHNICAL DATA

MINIMUM CONTINUOUS FLOW

REMARK: Due to manufacturing differences, the exhalation flow rate may differ from the values shown in the table ($\pm 10\%$).

Normative references	CE marking in accordance with European Directive 93/42/EEC, class IIa and subsequent amendments EN ISO 14971, ISO 5356-1, ISO 10993-1, EN ISO 17510-2
Resistance	at 50 lpm: 0,3 cmH ₂ O - at 100 lpm: 0,6 cmH ₂ O
Operating pressure	4-40 cmH ₂ O
Temperature range	5°C - 40°C
Dead space volume inside the mask	215 ml (Large)
Inspiratory and expiratory resistance during single fault condition when the anti-asphyxia valve is activated	0,6 cmH ₂ O
Anti-asphyxia valve opening pressure	0,3 cmH ₂ O
Anti-asphyxia valve closing pressure	0,1 cmH ₂ O
A-weighted sound pressure level according to ISO 3744	15,1 dB(A) \pm 1,5 dB(A)
A-weighted sound power level according to ISO 3744	23,1 dB(A) \pm 1,5 dB(A)

14 - LIMITED WARRANTY

This product by Air Liquide Medical Systems S.r.l. is guaranteed against material and manufacturing defects for a period of 120 days as from the date of purchase, provided it is used in accordance with the operating instructions and for the intended purpose only.

If the product fails under conditions of normal use, Air Liquide Medical Systems S.r.l. will repair or replace, at its option the

defective product or any of its components.

The manufacturer is liable for the safety, reliability and performance of the apparatus, provided that it is used in accordance with the instructions and for the intended purpose only, and that any repairs are made by the manufacturer or an approved service centre.

Throughout the warranty period keep the original packing.

Warranty claims on effective product must be made by the consumer contacting the dealer.

For further information on your warranty rights contact your dealer.

COPYRIGHT

The information contained in this booklet may not be used for anything other than the intended purpose. This booklet is the property of Air Liquide Medical Systems S.r.l. and no part of it may be reproduced without the manufacturer's written permission. All rights reserved.

TECHNICAL UPDATES

Air Liquide Medical Systems S.r.l. periodically reviews and modifies all its medical devices in order to improve their performance, safety and reliability. The instruction booklets are updated to ensure that they continue to be consistent with the characteristics of the devices marketed.

If the booklet accompanying the apparatus is mislaid, a replacement copy can be obtained from the Manufacturer by quoting the data shown on the label.

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AIR LIQUIDE MEDICAL SYSTEMS S.r.l.