Mode d'emploi / Instructions for use / Gebrauchsanweisung / Instrucciones de uso / Istruzioni per l'uso / Instruções de uso / Brugsanvisning / Instruksjoner for bruk / Användningsinstruktioner / Käyttöohjeet / Gebruiksaanwijzing / Instrukcja użycia / Návod k použití / Inštrukcie na používanie / Használati útmutató / Instruțiuni de folosire / Navodila za uporabo / Οδηγίες χρήσης / Инструкции по применению / Інструкція із застосування / 使用說明 / 使用說明書 / 사용지침 / تعليمات الاستخدام / הוראות שימוש



Respireo
Primo P

Masque avec coussin narinaire • Mask with nasal pads • Nasenpolster-Maske • Máscara con almohadillas nasales • Maschera a cuscinetti nasali • Máscara com almofadas nasais • Næsemaske med sidepudder • Neseputemaske • Mask med näskuddar • Nenätyynymskäsi • Neuskussenmasker • Maska z poduszczkami ochronnymi na nos • Maska s nosními polštářky • Maska s nosnými vankúšikmi • Orrpárnás maszk • Mască cu pernite nazale • Maska z nosnimi blazinicami • Μάσκα με μαξιλαράκι μύτης • Маска с назальными подушечками • Маска з назальними подушечками • 帶护鼻面罩 / 鼻パットマスク • 코 쿠션 마스크 • قناع بوساند أنفية • מסכה עם נחירונים


A


B


1


2


3


4

C


②


D

①

②

FRANÇAIS	p. 6	ROMÂNĂ.....	p. 47
ENGLISH.....	p. 8	SLOVENŠČINA.....	p. 50
DEUTSCH	p. 11	ΕΛΛΗΝΙΚΑ	p. 53
ESPAÑOL	p. 14	РУССКИЙ.....	p. 56
ITALIANO	p. 17	УКРАЇНСЬКА.....	p. 58
PORTUGUÊS	p. 20	繁體中文.....	p. 61
DANSK.....	p. 23	日本語.....	p. 63
NORSK	p. 25	한국어	p. 66
SVENSKA.....	p. 28	العربية.....	p. 70
SUOMI.....	p. 31	עברית.....	p. 73
NEDERLANDS.....	p. 34		
POLSKI	p. 36		
ČEŠTINA.....	p. 39		
SLOVENČINA.....	p. 42		
MAGYAR.....	p. 45		

1 - USAGE PRÉVU

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ENGLISH**1 - INTENDED USE**

This medical device belongs to the category "Masks for non-invasive ventilation and related connectors", that includes nasal, nasal pillow and oro-nasal masks, along with the related connectors, which are intended to be used at home or in a hospital setting by a single or by multiple patients who need for a non-invasive positive pressure ventilation (NIPPV) therapy.

INDICATION FOR USE AND INTENDED USERS

Respireo Primo P is a reusable nasal pillow mask with calibrated exhalation orifices (Vented), intended to be used at home or in a hospital setting by a single adult patient (weighing more than 30 kg) who has been prescribed for a non-invasive positive pressure ventilation therapy (NPPV), for example a continuous positive airway pressure therapy (CPAP) or Bilevel positive airway pressure therapy (BiPAP).

⚠ WARNINGS

- The mask should be used exclusively with positive pressure equipment recommended by your doctor or respiratory therapist supplying a minimum pressure of 4 cmH₂O: at lower pressures, exhaled air safety evacuation through the exhalation orifices is not guaranteed and so a partial re-inhalation may occur.
- Before using the mask, it needs to check that the exhalation orifices are not occluded, even partially, to not jeopardize the safety and the quality of the therapy.

- In case of discomfort, irritation or allergic reactions to any component of the mask, consult your doctor or respiratory physiotherapist.
- In case of additional oxygen administration, smoking or using open flames is prohibited.
- When using oxygen and the equipment is not in operation, turn off the oxygen dispenser.
- Do not use the mask in case of vomiting or nausea.
- Keep away from light.
- It is recommended to wash the mask before the first use.
- The mask shall be cleaned after each use.
- Frequency of cleaning, methods of cleaning or the use of cleaning agents other than those specified in this user manual can jeopardize the safety and the quality of the therapy.
- Before use, check the integrity of the mask. If damaged during transport, notify your dealer.
- In case of visible deterioration (cracks, tears, etc.), it needs to eliminate and replace the damaged mask component.
- Do not leave the components of the mask unattended, some of these may be swallowed by children.
- To ensure proper operation of the hose assembly, handle it with care, especially avoiding pulling it.
- Contact your dealer as a healthcare professional in case of functional problems, size problems or questions regarding the application of the medical device.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

2 - MASK COMPONENTS (section A)

- a) nasal pillow; b) frame; c) hose assembly; d) hose strap; e) rear buckle; f) top buckle.

REMARK: The pack contains three nasal pillows sizes S, M and L.

3 - WEARING THE MASK (section B)

- Insert the hose assembly (c) in the hose strap (d).
- Hold the frame (b) firmly in one hand, keeping the side pieces upwards. With the other hand, extend the stretch part of the frame with the other hand (B.1).
- Place the frame (b) over your head and insert the nasal pillow (a) into your nostrils (B.2).
- Adjust the tension of the rear strap (B.3).

- Adjust the tension of the top strap (B.4).
- Connect the swivel fitting of the hose assembly (c) to the equipment tube, then turn on the equipment at the prescribed pressure.

4 - USING THE MASK

To properly disassemble and reassemble the mask, refer to the respective illustrations in section C - D.

5 - CLEANING THE MASK

For daily cleaning, the mask must be disassembled in its components, each of which (except for the frame) must be carefully washed in warm water (about 30°C) using neutral soap and then rinsed thoroughly with clean fresh water. Finally, the parts must be air dried, away from direct sunlight.

To clean the frame, perform the same procedure weekly.

⚠ PRECAUTIONS

- Be very careful to remove any dirt from the exhalation holes.
- The frame can only be washed in water.
- Do not clean the mask with solutions containing alcohol, aromatic compounds, humidifiers, antibacterial agents, hydrating agents or aromatic oils as they would damage the product and reduce its life.
- Do not leave the mask immersed for more than 10 minutes.
- Do not clean the mask in the dishwasher.
- Check carefully all the components to make sure there are no soap deposits or foreign bodies inside the mask.

6 - TROUBLESHOOTING

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
The airflow doesn't reach the mask.	- Flow generator turned off or not working. - Breathing circuit ventilator not inserted correctly.	- Switch on or adjust the flow generator. - Connect the breathing circuit ventilator to the flow generator and the mask.

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
Noisy and annoying air leaks.	<ul style="list-style-type: none"> - The mask has not been assembled correctly. - The mask is not positioned correctly. - The size of the mask may not be correct. 	<ul style="list-style-type: none"> - Disassemble the mask and then reassemble it as described in the user manual. - Reposition the mask on your face and adjust the headgear straps as described in the user manual. - Contact your doctor or therapy technician to find the right size.
The exhalation orifices become dirty to the point of occlusion.	Incorrect or insufficient daily cleaning.	Clean the exhalation orifices using a soft-bristled toothbrush.


7 - DISPOSAL

The mask components are made of the materials indicated below. They do not contain hazardous substances and can therefore be disposed of with ordinary household waste.

COMPONENT	MATERIAL
Nasal pillow (a)	Silicone
Frame (b)	Nylon + Nylon/UBL/Polyurethane
Hose assembly (c)	Polycarbonate + Ethylene-Octene Copolymer + POM
Hose strap (d) and buckles (e), (f)	Nylon

8 - TECHNICAL DATA

PRESSURE - FLOW CURVE*

*Due to machining differences, the outflow from the exhalation orifices may vary ($\pm 10\%$).

Regulatory information	CE marking in accordance with Regulation 2017/745 (EU) and subsequent implementation amendments. Class IIa Main product standard UNI EN ISO 17510
Operating pressure range	4 cmH ₂ O / 24 cmH ₂ O
Operating temperature range	+5 °C / +40 °C
Dead space	Size S: 20 ml Size M: 21 ml Size L: 23 ml
Resistance	at 50 l/min: 0,6 cmH ₂ O at 100 l/min: 2,0 cmH ₂ O
A-Weighted Sound Pressure (UNI EN ISO 3744)	11,5 \pm 3,2 dB(A)
A-Weighted Sound Power (UNI EN ISO 3744)	18,2 \pm 3,2 dB(A)
Dimensions	40 x 55 x 70 mm
Weight	Size S: 65,5 g Size M: 66 g Size L: 66,5 g

9 - WARRANTY CERTIFICATE

Air Liquide Medical Systems S.r.l. guarantees its product against material or manufacturing defects for 180 days from the date of purchase, provided that the conditions of use in the instructions for use are observed.

If the product is faulty under normal use, Air Liquide Medical Systems S.r.l. will, at its discretion, repair or replace the faulty product or its components. The manufacturer shall be responsible for the safety, reliability and performance of the device if it is used in accordance with the operating instructions and its intended use.

Keep the original casing for the entire warranty period. The warranty in case of product defects can be applied by contacting your dealer.

For more information on the rights under the warranty, please contact your retailer.

COPYRIGHT

All the information in this manual may not be used for purposes other than the original. This manual is the property of Air Liquide Medical Systems S.r.l. and it may not be reproduced, in whole or in part, without written permission from the company. All rights are reserved.

TECHNICAL CHARACTERISTICS UPDATES

For a continuous improvement of performance, safety and reliability, all medical devices produced by Air Liquide Medical Systems S.r.l. are regularly subjected to review and change.

Instruction manuals are therefore modified to ensure their constant consistency with the characteristics of devices placed on the market.

If the instruction manual accompanying this device is lost, contact the manufacturer to get a copy of the version corresponding to the device by quoting the references shown on the label.

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AIR LIQUIDE MEDICAL SYSTEMS S.r.l.