MODE D'EMPLOI / INSTRUCTIONS FOR USE / GEBRAUCHSANWEISUNG / INSTRUCCIONES DE USO / ISTRUZIONI PER L'USO / INSTRUÇÕES DE USO  
BRUGSANVISNING / INSTRUKSJONER FOR BRUK / ANVÄNDNINGSSINSTRUKTIONER / KÄYTTÖOHJEET / GEBRUIKSAANWIJZING / INSTRUKCJA UŻYCIA  
NÁVOD K POUŽITÍ / INŠTRUKCIE NA POUŽÍVANIE / HASZNÁLATI ÚTMUTATÓ / INSTRUCTIUNI DE FOLOSIRE / NAVODILA ZA UPORABO / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ  
ИНСТРУКЦИИ ПО ПРИМЕНЕНИЮ / ІНСТРУКЦІЯ ІЗ ЗАСТОСУВАННЯ / 使用説明 / 使用説明書 / 사용 지침

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**1 - INTENDED USE**

This medical device belongs to the category "Masks for non-invasive ventilation and related connectors", that includes nasal, nasal pillow and oro-nasal masks, along with the related connectors, which are intended to be used at home or in a hospital setting by a single or by multiple patients who need for a non-invasive positive pressure ventilation (NIPPV) therapy.

**INDICATION FOR USE AND INTENDED USERS**

Respireo Hospital F Non Vented Disposable with anti-asphyxia valve is an oronasal mask without calibrated exhalation orifices and with an anti-suffocation valve, intended to be used for short time periods (less than 7 days) in a hospital setting by a single adult patient (weighing over 30kg) in controlled non-invasive positive pressure ventilation (NIPPV).

**⚠ WARNINGS**

- The mask should only be used with a ventilator fitted with an active exhalation valve and/or with a respiratory circuit with an exhalation valve recommended by your doctor or respiratory therapist.
- The anti-asphyxia valve fitting ensures patient safety even in case of failure of the positive pressure equipment. Do not use the mask if the valve is damaged or distorted.
- In case of discomfort, irritation or allergic reactions to any component of the mask, consult your doctor or respiratory physiotherapist
- In case of additional oxygen administration, smoking or using open flames is prohibited.
- When using oxygen and the equipment is not in operation, turn off the oxygen dispenser.
- Do not use the mask in case of vomiting or nausea.
- Keep away from light.
- Before use, check the integrity of the mask. If damaged during transport, notify your dealer.
- Do not leave the components of the mask unattended, some of these may be swallowed by children.
- The device is a single use device: the reuse on multiple patients may result in transmission of infections, serious harm or even the death of the patient.
- Contact your dealer as a healthcare professional in case of functional problems, size problems or questions regarding the

application of the medical device.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**2 - MASK COMPONENTS (section A)**

a) cushion; b) frame; c) forehead support; d) adjustment knob; e) sampling port cap; f) anti-asphyxia valve fitting with yellow connector; g) headgear; h) clip.

**3 - WEARING THE MASK (section B)**

- Turn the adjustment knob (d) until the forehead support (c) is fully extended (B.1).
- Holding the mask gently in close contact to the face, wear the headgear (g) over your head making sure that the cushion (a) adheres to the oronasal profile (B.2).
- Insert the two clips (h) into the respective seats of the frame (b) so that the lower straps pass under the ears (B.3).
- Slightly pull the straps to ensure the stability of the mask according to the sequence B.4.
- Connect the anti-asphyxia valve fitting with yellow connector (f) to the equipment tube (B.5), then turn on the equipment at the prescribed pressure.
- Lie down and turn the adjustment knob (d) to obtain a better adhesion to the face or better seal (B.6).
- If necessary, re-adjust the headgear.

**4 - REMOVING THE MASK (section C)**

To remove the mask, press the button on the front to release one of the clips (h). Lift the headgear over your head and remove the mask. The other straps do not need to be loosened nor disconnected.

**5 - TROUBLESHOOTING**

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
The airflow doesn't reach the mask.	<ul style="list-style-type: none"> <li>• Flow generator turned off or not working.</li> <li>• Breathing circuit ventilator not inserted correctly.</li> </ul>	<ul style="list-style-type: none"> <li>• Switch on or adjust the flow generator.</li> <li>• Connect the breathing circuit ventilator to the flow generator and the mask.</li> </ul>
Noisy and annoying air leaks.	<ul style="list-style-type: none"> <li>• The mask has not been assembled correctly.</li> <li>• The mask is not positioned correctly.</li> <li>• The size of the mask may not be correct.</li> </ul>	<ul style="list-style-type: none"> <li>• Disassemble the mask and then reassemble it as described in the user manual.</li> <li>• Reposition the mask on your face and adjust the headgear straps as described in the user manual.</li> <li>• Contact your doctor or therapy technician to find the right size.</li> </ul>

**6 - DISPOSAL**

The mask components are made of the materials indicated below. They do not contain hazardous substances and can therefore be disposed of with ordinary household waste.

COMPONENT	MATERIAL
Cushion (a)	Silicone
Frame (b)	Polycarbonate
Forehead support (c)	Silicone
Adjustment knob (d)	Polypropylene
Sampling port cap (e)	Silicone

COMPONENT	MATERIAL
Anti-asphyxia valve fitting with yellow connector (f)	Polycarbonate + Silicone
Headgear (g)	Nylon/UBL/Polyurethane
Clip (h)	Nylon

## 7 - TECHNICAL DATA

<b>Regulatory information</b>	CE marking in accordance with Regulation 2017/745 (EU) and subsequent implementation amendments. Class IIa Main product standard UNI EN ISO 5356-1
<b>Operating pressure range</b>	4 cmH <sub>2</sub> O / 32 cmH <sub>2</sub> O
<b>Operating temperature range</b>	+5 °C / +40 °C
<b>Opening pressure of the anti-asphyxia valve</b>	0,3 cmH <sub>2</sub> O
<b>Closing pressure of the anti-asphyxia valve</b>	0,1 cmH <sub>2</sub> O
<b>Dead space</b>	Size S: 283 mL Size M: 293 mL Size L: 329 mL
<b>Dimensions (without headgear)</b>	Size S: 220 x 100 x 120 mm Size M: 220 x 100 x 120 mm Size L: 230 x 100 x 120 mm
<b>Weight</b>	Size S: 163 g Size M: 165 g Size L: 166 g

## 8 - WARRANTY CERTIFICATE

Air Liquide Medical Systems S.r.l. guarantees its product against material or manufacturing defects for 180 days from the date of purchase, provided that the conditions of use in the instructions for use are observed.

If the product is faulty under normal use, Air Liquide Medical Systems S.r.l. will, at its discretion, repair or replace the faulty product or its components.

The manufacturer shall be responsible for the safety, reliability and performance of the device if it is used in accordance with the operating instructions and its intended use.

Keep the original casing for the entire warranty period. The warranty in case of product defects can be applied by contacting your dealer.

For more information on the rights under the warranty, please contact your retailer.

### COPYRIGHT

All the information in this manual may not be used for purposes other than the original. This manual is the property of Air Liquide Medical Systems S.r.l. and it may not be reproduced, in whole or in part, without written permission from the company. All rights are reserved.

### TECHNICAL CHARACTERISTICS UPDATES

For a continuous improvement of performance, safety and reliability, all medical devices produced by Air Liquide Medical Systems S.r.l. are regularly subjected to review and change.

Instruction manuals are therefore modified to ensure their constant consistency with the characteristics of devices placed on the market. If the instruction manual accompanying this device is lost, contact the manufacturer to get a copy of the version corresponding to the device by quoting the references shown on the label.

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AIR LIQUIDE MEDICAL SYSTEMS S.r.l.