Mode d'emploi / Instructions for use / Gebrauchsanweisung / Instrucciones de uso / Istruzioni per l'uso / Instruções de uso / Användningsinstruktioner /
Gebruiksaanwijzing / Instrukcja użycia / Návod k použití / Οδηγίες χρήσης / Инструкции по применению / 使用説明 / 使用説明書 / 사용 지침 / تعليمات الاستخدام


Respireo
Primo **F**




Masque oro-nasale avec fuites d'expiration / Oronasal mask with exhalation orifices / Mund-nasen-maske mit Ausatemungsöffnungen /
Mascarilla oro-nasal con orificios de exhalación / Maschera oronasale con fori di esalazione / Máscara oro-nasal com orificios de exalação
/ Oronasal mask med utandningshål / Oronasal masker met uitademingsopeningen / Maska oronasalna z otworami wentylacyjnymi
/ Oronasální maska s výdechovými otvory / Στοματορινική μάσκα με οπές εκπνοής / Оронозальная маска с вентиляционными
отверстиями / 带有标准通气口鼻面罩 / 鼻と口のマスク通気穴付き / 호기 구멍이 있는 페이스얼 마스크 / قناع فموي انفي مع فتحات للزفير

A



B


B



④

⑤

⑥

C**D**

FRANÇAIS	p. 09
ENGLISH.....	p. 12
DEUTSCH	p. 15
ESPAÑOL	p. 19
ITALIANO	p. 23
PORTUGUÊS	p. 26
SVENSKA.....	p. 29
NEDERLANDS.....	p. 33
POLSKI.....	p. 37
ČEŠTINA.....	p. 41
ΕΛΛΗΝΙΚΑ	p. 44
РУССКИЙ.....	p. 48
繁體中文.....	p. 52
日本語.....	p. 55
한국어.....	p. 58
العربية.....	p. 64

1 - USAGE PRÉVU

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ENGLISH

1 - INTENDED USE

This medical device belongs to the category "Masks for non-invasive ventilation and related connectors", that includes nasal, nasal pillow and oro-nasal masks, along with the related connectors, which are intended to be used at home or in a hospital setting by a single or by multiple patients who need for a non-invasive positive pressure ventilation (NIPPV) therapy.

INDICATION FOR USE AND INTENDED USERS

Respireo Primo F is a reusable oronasal mask with calibrated exhalation orifices (Vented), intended to be used at home or in a hospital setting by a

ENGLISH

single adult patient (weighing more than 30 kg) who has been prescribed for a non-invasive positive pressure ventilation therapy (NPPV), for example a continuous positive airway pressure therapy (CPAP) or Bilevel positive airway pressure therapy (BiPAP).

WARNINGS

- The mask should be used exclusively with positive pressure equipment recommended by your doctor or respiratory therapist supplying a minimum pressure of 4 cmH₂O: at lower pressures, exhaled air safety evacuation through the exhalation orifices is not guaranteed and so a partial re-inhalation may occur.
- The anti-asphyxia valve fitting ensures patient safety even in case of failure of the positive pressure equipment. Do not use the mask if the valve is damaged or distorted.
- Before using the mask, it needs to check that the exhalation orifices are not occluded, even partially, to not jeopardize the safety and the quality of the therapy.
- In case of discomfort, irritation or allergic reactions to any component of the mask, consult your doctor or respiratory physiotherapist.
- In case of additional oxygen administration, smoking or using open flames is prohibited.
- When using oxygen and the equipment is not in operation, turn off the oxygen dispenser.
- Do not use the mask in case of vomiting or nausea.
- Keep away from light.
- It is recommended to wash the mask before the first use.
- The mask shall be cleaned after each use.
- Frequency of cleaning, methods of cleaning or the use of cleaning agents other than those specified in this user manual, can jeopardize the safety and the quality of the therapy.
- Before use, check the integrity of the mask. If damaged during transport, notify your dealer.
- In case of visible deterioration (cracks, tears, etc.), it needs to eliminate and replace the damaged mask component.
- Do not leave the components of the mask unattended, some of these may be swallowed by children.

- Contact your dealer as a healthcare professional in case of functional problems, size problems or questions regarding the application of the medical device.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

2 - MASK COMPONENTS (section A)

a) cushion; b) frame; c) forehead support; d) adjustment knob; e) sampling port cap; f) anti-asphyxia valve fitting; g) swivel fitting; h) headgear. i) clip.

3 - WEARING THE MASK (section B)

- Turn the adjustment knob (d) until the forehead support is fully extended (B.1).
- Holding the mask gently in close contact to the face, wear the headgear (h) over your head making sure that the cushion (a) adheres to the oronasal profile (B.2).
- Insert the two clips (i) into the respective seats of the frame (b) so that the lower straps pass under the ears (B.3).
- Slightly pull the straps to ensure the stability of the mask according to the sequence B.4.
- Connect the swivel fitting (g) to the equipment tube (B.5), then turn on the equipment at the prescribed pressure.
- Lie down and turn the adjustment knob (d) to obtain a better adhesion to the face or better seal (B.6).
- If necessary, re-adjust the headgear.

4 - REMOVING THE MASK (section C)

To remove the mask, press the button on the front to release one of the clips (i). Lift the headgear over your head and remove the mask. The other straps do not need to be loosened nor disconnected.

5 - DISASSEMBLING THE MASK

- Disconnect the breathing circuit ventilator from the swivel fitting (g).
- Remove the two clips (i) from the frame (b).
- Remove the headgear (h) from the frame (b) and the clips (i) by opening and removing the four straps.
- Remove the cushion (a), the forehead support (c), the port cap (e) and the blue elbow connector (f) from the frame (b).

- Detach the swivel fitting (g) from the anti-asphyxia valve fitting (f).

6 - REASSEMBLING THE MASK

- Connect the swivel fitting (g) to the anti-asphyxia valve fitting (f).
- Connect the anti-asphyxia valve fitting (f) into the frame (b) by applying light pressure.
- Insert the forehead support (c) into the appropriate seat of the frame (b) keeping the curved edge upwards.
- Insert the grooved side of the cushion (a) into the frame (b) starting from the upper top of the cushion and continuing along the contour of the mask pushing until the stop segments positioned on the frame are reached (section D).
- Insert the port cap (e) into its seat.
- Insert the two upper headgear (h) straps into the respective frame (b) slots and both the lower straps into the clip (i) slots.
- Insert the two clips (i) into the respective seats of the frame (b).

7 - CLEANING THE MASK

For daily cleaning, the mask must be disassembled in its components, each of which (except for the headgear) must be carefully washed in warm water (about 30°C) using neutral soap and then rinsed thoroughly with clean fresh water. Finally, the parts must be air dried, away from direct sunlight.

To clean the headgear, perform the same procedure weekly.

For a deeper decontamination of the mask, consult the “Guide to decontamination, disinfection and sterilization” available electronically on the website it.medicaldevice.airliquide.com or at your dealer, together with the instructions provided by the manufacturer of the disinfectant solution, especially regarding the compatibility with the materials of the various components.

⚠ PRECAUTIONS

- Note the position of the adjustment knob at the first use to facilitate the ideal repositioning of the mask at each subsequent use.
- Remove the mask after disconnecting at least one of the lower headgear straps.
- Do not open the straps to avoid having to adjust the headgear again the next time you use it.
- Be very careful to remove any dirt from the exhalation holes.

- The headgear can only be washed in water.
- Do not clean the mask with solutions containing alcohol, aromatic compounds, humidifiers, antibacterial agents, hydrating agents or aromatic oils as they would damage the product and reduce its life.
- Do not leave the mask immersed for more than 10 minutes.
- Do not clean the mask in the dishwasher.
- Check carefully all the components to make sure there are no soap deposits or foreign bodies inside the mask.

8 - TROUBLESHOOTING

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
The airflow doesn't reach the mask.	- Flow generator turned off or not working. - Breathing circuit ventilator not inserted correctly.	- Switch on or adjust the flow generator. - Connect the breathing circuit ventilator to the flow generator and the mask.
Noisy and annoying air leaks.	- The mask has not been assembled correctly. - The mask is not positioned correctly. - The size of the mask may not be correct.	- Disassemble the mask and then reassemble it as described in the user manual. - Reposition the mask on your face and adjust the headgear straps as described in the user manual. - Contact your doctor or therapy technician to find the right size.
The exhalation orifices become dirty to the point of occlusion.	Incorrect or insufficient daily cleaning.	Clean the exhalation orifices using a soft-bristled toothbrush.


9 - DISPOSAL

The mask components are made of the materials indicated below. They do not contain hazardous substances and can therefore be disposed of with ordinary household waste.

COMPONENT	MATERIAL
Cushion (a)	Silicone
Frame (b)	Polycarbonate
Forehead support (c)	Silicone
Adjustment knob (d)	Polypropylene
Sampling port cap (e)	Silicone
Anti-asphyxia valve fitting (f)	Polycarbonate + Silicone
Swivel fitting (g)	Polycarbonate
Headgear (h)	Nylon/UBL/Polyurethane
Clip (i)	Nylon

10 - TECHNICAL DATA

PRESSURE - FLOW CURVE*

*Due to machining differences, the outflow from the exhalation orifices may vary ($\pm 10\%$).

Regulatory information	CE marking in accordance with Regulation 2017/745 (EU) and subsequent implementation amendments. Class IIa Main product standard UNI EN ISO 17510
Operating pressure range	4 cmH ₂ O / 32 cmH ₂ O
Operating temperature range	+5 °C / +40 °C

ENGLISH

Dead space	Size S: 285 mL Size M: 309 mL Size L: 334 mL
Resistance	at 50 l/min: 0,3 cmH ₂ O at 100 l/min: 0,6 cmH ₂ O
Resistance in single fault condition, when the anti-asphyxia valve is opened in the atmosphere	Inhalation at 50 l/min: 0,6 cmH ₂ O Exhalation at 50 l/min: 0,6 cmH ₂ O
Opening pressure of the anti-asphyxia valve	0,3 cmH ₂ O
Closing pressure of the anti-asphyxia valve	0,1 cmH ₂ O
A-Weighted Sound Pressure (UNI EN ISO 3744)	15,1 ± 1,5 dB(A)
A-Weighted Sound Power (UNI EN ISO 3744)	23,1 ± 1,5 dB(A)
Dimensions (without headgear)	Size S: 220 x 100 x 120 mm Size M: 220 x 100 x 120 mm Size L: 230 x 100 x 120 mm
Weight	Size S: 175 g Size M: 175,5 g Size L: 176 g

11 - WARRANTY CERTIFICATE

Air Liquide Medical Systems S.r.l. guarantees its product against material or manufacturing defects for 180 days from the date of purchase, provided that the conditions of use in the instructions for use are observed.

If the product is faulty under normal use, Air Liquide Medical Systems S.r.l. will, at its discretion, repair or replace the faulty product or its components.

The manufacturer shall be responsible for the safety, reliability and performance of the device if it is used in accordance with the operating instructions and its intended use.

Keep the original casing for the entire warranty period. The warranty in case of product defects can be applied by contacting your dealer.

For more information on the rights under the warranty, please contact your retailer.

COPYRIGHT

All the information in this manual may not be used for purposes other than the original. This manual is the property of Air Liquide Medical Systems S.r.l.

ENGLISH

and it may not be reproduced, in whole or in part, without written permission from the company. All rights are reserved.

TECHNICAL CHARACTERISTICS UPDATES

For a continuous improvement of performance, safety and reliability, all medical devices produced by Air Liquide Medical Systems S.r.l. are regularly subjected to review and change.

Instruction manuals are therefore modified to ensure their constant consistency with the characteristics of devices placed on the market.

If the instruction manual accompanying this device is lost, contact the manufacturer to get a copy of the version corresponding to the device by quoting the references shown on the label.

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AIR LIQUIDE MEDICAL SYSTEMS S.r.l.