Mode d'emploi / Instructions for use / Gebrauchsanweisung / Instrucciones de uso / Istruzioni per l'uso / Instruções de uso / Användningsinstruktioner /
Gebruiksaanwijzing / Instrukcja użycia / Návod k použití / Οδηγίες χρήσης / Инструкция по использованию / 使用说明 / 使用説明書 / 사용 지침 / تعليمات الاستخدام

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



Masque Oro-nasale / Oronasal Mask / Mund-Nasen-Maske / Mascarilla Oronasal / Maschera Oronasale / Máscara Oro-Nasal /
Oronasal Mask / Mondneusmasker / Maska ustno-nosowa / Oronazální maska / Στοματορινική Μάσκα / Оronазальная Маска / 口鼻罩
/ 口鼻マスク / 구비강 마스크 / قناع فموي أنفي






A

B**1****2****3**

B

C**D**

FRANÇAIS	7
ENGLISH.....	10
DEUTSCH	12
ESPAÑOL.....	16
ITALIANO	19
PORTUGUÊS.....	21
SVENSKA.....	25
NEDERLANDS	27
POLSKI.....	31
ČEŠTINA	34
РУССКИЙ	40
简体中文	43
日本語	46
한국어.....	48
ةيبرعلا	53

1 - USAGE PRÉVU

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1 - INTENDED USE

This medical device belongs to the category "Masks for non-invasive ventilation and related connectors", that includes nasal, nasal pillow and oro-nasal masks, along with the related connectors, which are intended to be used at home or in a hospital setting by a single or by multiple patients who need for a non-invasive positive pressure ventilation (NIPPV) therapy.

INDICATION FOR USE AND INTENDED USERS

Respireo Primo F Non Vented is a reusable oronasal mask without calibrated exhalation holes intended to be used at home or in a hospital setting by a single adult patient (weighing more than 30 kg) under Non-invasive Positive-Pressure Ventilation (NPPV).

⚠ WARNINGS

- The mask should only be used with a ventilator fitted with an active exhalation valve and/or with a respiratory circuit with an exhalation valve recommended by your doctor or respiratory therapist.
- In case of discomfort, irritation or allergic reactions to any component of the mask, consult your doctor or respiratory physiotherapist
- In case of additional oxygen administration, smoking or using open flames is prohibited.
- When using oxygen and the equipment is not in operation, turn off the oxygen dispenser.
- Do not use the mask in case of vomiting or nausea.
- Keep away from light.
- It is recommended to wash the mask before the first use.
- The mask shall be cleaned after each use.
- Frequency of cleaning, methods of cleaning or the use of cleaning agents other than those specified in this user manual, can jeopardize the safety and the quality of the therapy.
- Before use, check the integrity of the mask. If damaged during transport, notify your dealer.
- In case of visible deterioration (cracks, tears, etc.), it needs to eliminate and replace the damaged mask component.
- Do not leave the components of the mask unattended, some of these may be swallowed by children.

- Contact your dealer as a healthcare professional in case of functional problems, size problems or questions regarding the application of the medical device.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

2 - MASK COMPONENTS (section A)

a) cushion; b) frame; c) forehead support; d) adjustment knob; e) sampling port cap; f) blue elbow connector; g) headgear; h) clip.

3 - WEARING THE MASK (section B)

- Turn the adjustment knob (d) until the forehead support (c) is fully extended (B.1).
- Holding the mask gently in close contact to the face, wear the headgear (g) over your head making sure that the cushion (a) adheres to the oronasal profile (B.2).
- Insert the two clips (h) into the respective seats of the frame (b) so that the lower straps pass under the ears (B.3).
- Slightly pull the straps to ensure the stability of the mask according to the sequence B.4.
- Connect the blue elbow connector (f) to the equipment tube (B.5), then turn on the equipment at the prescribed pressure.
- Lie down and turn the adjustment knob (d) to obtain a better adhesion to the face or better seal (B.6).
- If necessary, re-adjust the headgear.

4 - REMOVING THE MASK (section C)

To remove the mask, press the button on the front to release one of the clips (h). Lift the headgear over your head and remove the mask. The other straps do not need to be loosened nor disconnected.

5 - DISASSEMBLING THE MASK

- Disconnect the breathing circuit ventilator from the blue elbow connector (f).
- Remove the two clips (h) from the frame (b).
- Remove the headgear (g) from the frame (b) and the clips (h) by opening and removing the four straps.

- Remove the cushion (a), the forehead support (c), the port cap (e) and the blue elbow connector (f) from the frame (b).

6 - REASSEMBLING THE MASK

- Connect the blue elbow connector (f) into the frame (b) by applying light pressure.
- Insert the forehead support (c) into the appropriate seat of the frame (b) keeping the curved edge upwards.
- Insert the grooved side of the cushion (a) into the frame (b) starting from the upper top of the cushion and continuing along the contour of the mask pushing until the stop segments positioned on the frame are reached (section D).
- Insert the port cap (e) into its seat.
- Insert the two upper headgear (g) straps into the respective frame (b) slots and both the lower straps into the clip (h) slots.
- Insert the two clips (h) into the respective seats of the frame (b).

7 - CLEANING THE MASK

For daily cleaning, the mask must be disassembled in its components, each of which (except for the headgear) must be carefully washed in warm water (about 30°C) using neutral soap and then rinsed thoroughly with clean fresh water. Finally, the parts must be air dried, away from direct sunlight.

To clean the headgear, perform the same procedure weekly.

For a deeper decontamination of the mask, consult the “Guide to decontamination, disinfection and sterilization” available electronically on the website medicaldevice.airliquide.com or at your dealer, together with the instructions provided by the manufacturer of the disinfectant solution, especially regarding the compatibility with the materials of the various components.

⚠ PRECAUTIONS

- Note the position of the adjustment knob at the first use to facilitate the ideal repositioning of the mask at each subsequent use.
- The headgear can only be washed in water.
- Do not clean the mask with solutions containing alcohol, aromatic compounds, humidifiers, antibacterial agents, hydrating agents or aromatic oils as they would damage the product and reduce its life.

- Do not leave the mask immersed for more than 10 minutes.
- Do not clean the mask in the dishwasher.
- Check carefully all the components to make sure there are no soap deposits or foreign bodies inside the mask.

8 - TROUBLESHOOTING

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
The airflow doesn't reach the mask.	<ul style="list-style-type: none"> - Flow generator turned off or not working. - Breathing circuit ventilator not inserted correctly. 	<ul style="list-style-type: none"> - Switch on or adjust the flow generator. - Connect the breathing circuit ventilator to the flow generator and the mask.
Noisy and annoying air leaks.	<ul style="list-style-type: none"> - The mask has not been assembled correctly. - The mask is not positioned correctly. - The size of the mask may not be correct. 	<ul style="list-style-type: none"> - Disassemble the mask and then reassemble it as described in the user manual. - Reposition the mask on your face and adjust the headgear straps as described in the user manual. - Contact your doctor or therapy technician to find the right size.

9 - DISPOSAL

The mask components are made of the materials indicated below. They do not contain hazardous substances and can therefore be disposed of with ordinary household waste.

COMPONENT	MATERIAL
Cushion (a)	Silicone
Frame (b)	Polycarbonate
Forehead support (c)	Silicone

ENGLISH

COMPONENT	MATERIAL
Adjustment knob (d)	Polypropylene
Sampling port cap (e)	Silicone
Blue elbow connector (f)	Polycarbonate
Headgear (g)	Nylon/UBL/Polyurethane
Clip (h)	Nylon

10 - TECHNICAL DATA

Regulatory information	CE marking in accordance with Regulation 2017/745 (EU) and subsequent implementation amendments. Class IIa Main product standard UNI EN ISO 5356-1
Operating pressure range	4 cmH ₂ O / 32 cmH ₂ O
Operating temperature range	+5°C / +40°C
Dead space	Size S: 283 mL Size M: 308 mL Size L: 329 mL
Dimensions (without headgear)	Size S: 205 x 100 x 120 mm Size M: 205 x 100 x 120 mm Size L: 210 x 100 x 120 mm
Weight	Size S: 168.5 g Size M: 169 g Size L: 169.5 g

11 - WARRANTY CERTIFICATE

Air Liquide Medical Systems S.r.l. guarantees its product against material or manufacturing defects for 180 days from the date of purchase, provided that the conditions of use in the instructions for use are observed. If the product is faulty under normal use, Air Liquide Medical Systems S.r.l. will, at its discretion, repair or replace the faulty product or its components.

ENGLISH

The manufacturer shall be responsible for the safety, reliability and performance of the device if it is used in accordance with the operating instructions and its intended use. Keep the original casing for the entire warranty period. The warranty in case of product defects can be applied by quoting your dealer. For more information on the rights under the warranty, please contact your retailer.

COPYRIGHT

All the information in this manual may not be used for purposes other than the original. This manual is the property of Air Liquide Medical Systems S.r.l. and it may not be reproduced, in whole or in part, without written permission from the company. All rights are reserved.

TECHNICAL CHARACTERISTICS UPDATES

For a continuous improvement of performance, safety and reliability, all medical devices produced by Air Liquide Medical Systems S.r.l. are regularly subjected to review and change.

Instruction manuals are therefore modified to ensure their constant consistency with the characteristics of devices placed on the market. If the instruction manual accompanying this device is lost, contact the manufacturer to get a copy of the version corresponding to the device by quoting the references shown on the label.manufacturer by quoting the data shown on the rating plate.

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AIR LIQUIDE MEDICAL SYSTEMS S.r.l.