MODE D'EMPLOI / INSTRUCTIONS FOR USE / GEBRAUCHSANWEISUNG / INSTRUCCIONES DE USO / ISTRUZIONI PER L'USO / INSTRUÇÕES DE USO /
BRUGSANVISNING / INSTRUKSJONER FOR BRUK / ANVÄNDNINGSSINSTRUKTIONER / KÄYTTÖOHJEET / GEBRUIKSAANWIJZING / INSTRUKCJA UŻYCIA /
NÁVOD K POUŽITÍ / INŠTRUKCIE NA POUŽÍVANIE / HASZNÁLATI ÚTMUTATÓ / INSTRUCCIUNI DE FOLOSIRE / NAVODILA ZA UPORABO / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
/ ИНСТРУКЦИИ ПО ПРИМЕНЕНИЮ / ІНСТРУКЦІЯ ІЗ ЗАСТОСУВАННЯ / 使用説明 / 使用説明書 / 사용 지침 / تعليمات الاستخدام

Respireo


Hospital F

Vented Disposable


Masque oro-nasale avec fuites d'expiration / Oronasal mask with exhalation orifices / Mund-nasen-maske mit Ausatemungsöffnungen /
Mascarilla oro-nasal con orificios de exhalación / Maschera oronasale con fori di esalazione / Máscara oro-nasal com orificios de exalação
/ Oronasal maske med udåndingsåbninger / Oronasal maske som er uten utpustningsåpninger / Oronasal mask med utandningshål /
Oronasaalinen kasvomaski jossa on uloshengityksen aukot / Oronasaal masker met uitademingsopening / Maska oronasalna z otworami
wentylacyjnymi / Oronasální maska s výdechovými otvory / Oronazálna maska s výdychovými otvormi / Oronasal maszk kilégzőnyílásokkal
/ Mască oronazală cu orificii de dispersare a aerului expirat / Oronazalna maska s luknjami za izdihavanje / Στοματορινική μάσκα με οπές
εκπνοής / Ороназальная маска с вентиляционными отверстиями / Ороназальна маска з вентиляційними отворами / 带有标准通气
口鼻面罩 / 鼻と口のマスク通気穴付き / 호기 구멍이 있는 페이스 마스크 / قناع فموي أنفي مع فتحات للزفير

A


B


1


2





3

4

B**5****6****C**

FRANÇAIS	7	ČEŠTINA	42
ENGLISH	10	SLOVENČINA	45
DEUTSCH	13	MAGYAR	48
ESPAÑOL	16	ROMÂNĂ	51
ITALIANO	19	SLOVENŠČINA	54
PORTUGUÊS	22	ΕΛΛΗΝΙΚΑ	57
DANSK	25	РУССКИЙ	60
NORSK	28	УКРАЇНЕЦЬ	64
SVENSKA	31	繁體中文	67
SUOMI	34	日本語	69
NEDERLANDS	36	한국어	71
POLSKI	39	العربية	76

1 - USAGE PRÉVU

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1 - INTENDED USE

This medical device belongs to the category “Masks for non-invasive ventilation and related connectors”, that includes nasal, nasal pillow and oro-nasal masks, along with the related connectors, which are intended to be used at home or in a hospital setting by a single or by multiple patients who need for a non-invasive positive pressure ventilation (NIPPV) therapy.

INDICATION FOR USE AND INTENDED USERS

Respireo Hospital F Vented Disposable is an oronasal mask with calibrated exhalation orifices intended to be used for short time periods (less than 7 days) in a hospital setting by a single adult patient (weighing over 30kg), who has been prescribed for a non-invasive positive pressure ventilation therapy (NIPPV), e.g. continuous positive airway pressure therapy (CPAP) or Bilevel positive airway pressure therapy (BiPAP).

⚠ WARNINGS

- The mask should be used exclusively with positive pressure equipment recommended by your doctor or respiratory therapist supplying a minimum pressure of 4 cmH₂O: at lower pressures, exhaled air safety evacuation through the exhalation orifices is not guaranteed and so a partial re-inhalation may occur.
- The anti-asphyxia valve fitting ensures patient safety even in case of failure of the positive pressure equipment. Do not use the mask if the valve is damaged or distorted.
- Before using the mask, it needs to check that the exhalation orifices are not occluded, even partially, to not jeopardize the safety and the quality of the therapy.
- In case of discomfort, irritation or allergic reactions to any component of the mask, consult your doctor or respiratory physiotherapist
- In case of additional oxygen administration, smoking or using open flames is prohibited.
- When using oxygen and the equipment is not in operation, turn off the oxygen dispenser.
- Do not use the mask in case of vomiting or nausea.
- Keep away from light.

- Before use, check the integrity of the mask. If damaged during transport, notify your dealer.
- Do not leave the components of the mask unattended, some of these may be swallowed by children.
- The device is a single use device: the reuse on multiple patients may result in transmission of infections, serious harm or even the death of the patient.
- Contact your dealer as a healthcare professional in case of functional problems, size problems or questions regarding the application of the medical device.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

2 - MASK COMPONENTS (section A)

a) cushion; b) frame; c) forehead support; d) adjustment knob; e) sampling port cap; f) anti-asphyxia valve fitting; g) swivel fitting; h) headgear; i) clip.

3 - WEARING THE MASK (section B)

- Turn the adjustment knob (d) until the forehead support (c) is fully extended (B.1).
- Holding the mask gently in close contact to the face, wear the headgear (h) over your head making sure that the cushion (a) adheres to the oronasal profile (B.2).
- Insert the two clips (i) into the respective seats of the frame (b) so that the lower straps pass under the ears (B.3).
- Slightly pull the straps to ensure the stability of the mask according to the sequence B.4.
- Connect the swivel fitting (g) to the equipment tube (B.5), then turn on the equipment at the prescribed pressure.
- Lie down and turn the adjustment knob (d) to obtain a better adhesion to the face or better seal (B.6).
- If necessary, re-adjust the headgear.

4 - REMOVING THE MASK (section C)

To remove the mask, press the button on the front to release one of the clips (i). Lift the headgear over your head and remove the mask. The other straps do not need to be loosened nor disconnected.

5 - TROUBLESHOOTING

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
The airflow doesn't reach the mask.	<ul style="list-style-type: none"> Flow generator turned off or not working. Breathing circuit ventilator not inserted correctly. 	<ul style="list-style-type: none"> Switch on or adjust the flow generator. Connect the breathing circuit ventilator to the flow generator and the mask.
Noisy and annoying air leaks.	<ul style="list-style-type: none"> The mask has not been assembled correctly. The mask is not positioned correctly. The size of the mask may not be correct. 	<ul style="list-style-type: none"> Disassemble the mask and then reassemble it as described in the user manual. Reposition the mask on your face and adjust the headgear straps as described in the user manual. Contact your doctor or therapy technician to find the right size.

6 - DISPOSAL


The mask components are made of the materials indicated below. They do not contain hazardous substances and can therefore be disposed of with ordinary household waste.

COMPONENT	MATERIAL
Cushion (a)	Silicone
Frame (b)	Polycarbonate

COMPONENT	MATERIAL
Forehead support (c)	Silicone
Adjustment knob (d)	Polypropylene
Sampling port cap (e)	Silicone
Anti-asphyxia valve fitting (f)	Polycarbonate + Silicone
Swivel fitting (g)	Polycarbonate
Headgear (h)	Nylon/UBL/Polyurethane
Clip (i)	Nylon

7 - TECHNICAL DATA

PRESSURE - FLOW CURVE*

*Due to machining differences, the outflow from the exhalation orifices may vary ($\pm 10\%$).

Regulatory information	
	CE marking in accordance with Regulation 2017/745 (EU) and subsequent implementation amendments. Class IIa Main product standard UNI EN ISO 17510

Operating pressure range	4 cmH ₂ O / 32 cmH ₂ O
Operating temperature range	+5 °C / +40 °C
Dead space	Size S: 285 mL Size M: 310 mL Size L: 332 mL
Resistance	at 50 l/min: 0,3 cmH ₂ O at 100 l/min: 0,6 cmH ₂ O
Resistance in single fault condition, when the anti-asphyxia valve is opened in the atmosphere	Inhalation at 50 l/min: 0,6 cmH ₂ O Exhalation at 50 l/min: 0,6 cmH ₂ O
Opening pressure of the anti-asphyxia valve	0,3 cmH ₂ O
Closing pressure of the anti-asphyxia valve	0,1 cmH ₂ O
A-Weighted Sound Pressure (UNI EN ISO 3744)	15,1 ± 3,2 dB(A)
A-Weighted Sound Power (UNI EN ISO 3744)	23,1 ± 3,2 dB(A)
Dimensions (without headgear)	Size S: 220 x 100 x 120 mm Size M: 220 x 100 x 120 mm Size L: 230 x 100 x 120 mm
Weight	Size S: 165 g Size M: 167 g Size L: 168 g

8 - WARRANTY CERTIFICATE

Air Liquide Medical Systems S.r.l. guarantees its product against material or manufacturing defects for 180 days from the date of purchase, provided that the conditions of use in the instructions for use are observed.

If the product is faulty under normal use, Air Liquide Medical Systems S.r.l. will, at its discretion, repair or replace the faulty product or its components.

The manufacturer shall be responsible for the safety, reliability and performance of the device if it is used in accordance with the operating instructions and its intended use.

Keep the original casing for the entire warranty period. The warranty in case of product defects can be applied by contacting your dealer.

For more information on the rights under the warranty, please contact your retailer.

COPYRIGHT

All the information in this manual may not be used for purposes other than the original. This manual is the property of Air Liquide Medical Systems S.r.l. and it may not be reproduced, in whole or in part, without written permission from the company. All rights are reserved.

TECHNICAL CHARACTERISTICS UPDATES

For a continuous improvement of performance, safety and reliability, all medical devices produced by Air Liquide Medical Systems S.r.l. are regularly subjected to review and change.

Instruction manuals are therefore modified to ensure their constant consistency with the characteristics of devices placed on the market. If the instruction manual accompanying this device is lost, contact the manufacturer to get a copy of the version corresponding to the device by quoting the references shown on the label.

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AIR LIQUIDE MEDICAL SYSTEMS S.r.l.