Mode d'emploi / Instructions for use / Gebrauchsanweisung / Instrucciones de uso / Istruzioni per l'uso / Instruções de uso / Brugsanvisning  
Instruksjoner for bruk / Användningsinstruktioner / Käyttöohjeet / Gebruiksaanwijzing / Instrukcja użycia / Návod k použití / Návod na použitie / Használati útmutató  
Instrucțiuni de utilizare / Navodila za uporabo / Οδηγίες χρήσης / Инструкция по использованию / Інструкція із застосування / 使用说明 / 使用說明書 / 사용 지침


# AInest N1 Silent

Masque nasal avec fuites d'expiration / Nasal mask with exhalation orifices / Nasenmaske mit Ausatemöffnungen / Mascarilla nasal con orificios de exhalación / Maschera nasale con fori di esalazione / Máscara nasal com orificios de exalação / Nasalmaske med udåndingsåbninger / Nesemaske med utpustningsåpninger / Näsmask med utandningshål / Nenämäski uloshengitysaukoilla / Neusmasker met uitademingsopeningen / Maska nosowa z otworami wydechowymi / Nosní maska s výdechovými otvory / Nazálna maska s výdychovými otvorami / Orrmaszk kilégzőnyílásokkal / Mască nazală cu orificii de expirație / Nosna maska z luknjami za izdihavanje / Ρινική μάσκα με οπές εκπνοής / Назальная маска с вентиляционными отверстиями / Назальна маска з вентиляційними отворами / 带有标准排气孔的鼻罩 / 呼吸排出孔付き鼻マスク / 호기구가 있는 비강 마스크


**A**


**B**

①

②


③


**B**


④







⑤





⑥

**C****1****2****D****1****3****2****4**

---

**E****1****2****3****4**

FRANÇAIS .....	8	ČEŠTINA .....	53
ENGLISH .....	11	SLOVENČINA .....	56
DEUTSCH .....	15	MAGYAR .....	60
ESPAÑOL .....	19	ROMÂNĂ .....	64
ITALIANO .....	23	SLOVENŠČINA .....	68
PORTUGUÊS .....	26	ΕΛΛΗΝΙΚΑ .....	71
DANSK .....	30	РУССКИЙ .....	75
NORSK .....	34	УКРАЇНСЬКА .....	79
SVENSKA .....	38	简体中文 .....	83
SUOMI .....	41	日本語 .....	86
NEDERLANDS .....	45	한국어 .....	89
POLSKI .....	49		

## 1 - USAGE PRÉVU

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## ENGLISH

### 1 - INTENDED USE

This medical device belongs to the category "Masks for non-invasive ventilation and related connectors", that includes nasal, nasal pillow and oro-nasal masks, along with the related connectors, which are intended to be used at home or in a hospital setting by a single or by multiple patients who need for a non-invasive positive pressure ventilation (NIPPV) therapy.

### INDICATION FOR USE AND INTENDED USERS

Alnest N1 Silent is a reusable nasal mask with calibrated exhalation orifices (Vented), intended to be used at home or in a hospital setting by multiple adult patients (weighing more than 30 kg) who have been prescribed for a non-invasive positive pressure ventilation therapy (NPPV), for example a continuous positive airway pressure therapy (CPAP) or Bilevel positive airway pressure therapy (BiPAP).

**⚠ WARNINGS**

- The mask should be used exclusively with positive pressure equipment recommended by your doctor or respiratory therapist supplying a minimum pressure of 4 cmH<sub>2</sub>O: at lower pressures, exhaled air safety evacuation through the exhalation orifices is not guaranteed and so a partial re-inhalation may occur.
- Before using the mask, it needs to check that the exhalation orifices are not occluded, even partially, to not jeopardize the safety and the quality of the therapy.
- In case of discomfort, irritation or allergic reactions to any component of the mask, consult your doctor or respiratory physiotherapist.
- In case of additional oxygen administration, smoking or using open flames is prohibited.
- When using oxygen and the equipment is not in operation, turn off the oxygen dispenser.
- Do not use the mask in case of vomiting or nausea.
- Keep away from light.
- It is recommended to wash the mask before the first use.
- The mask shall be cleaned after each use.
- Frequency of cleaning, methods of cleaning or the use of cleaning agents other than those specified in this user manual, or exceeding the recommended number of regeneration cycles, can jeopardize the safety and the quality of the therapy.
- Before use, check the integrity of the mask. If damaged during transport, notify your dealer.
- In case of visible deterioration (cracks, tears, etc.), it needs to eliminate and replace the damaged mask component.
- Do not leave the components of the mask unattended, some of these may be swallowed by children.
- Contact your dealer as a healthcare professional in case of functional problems, size problems or questions regarding the application of the medical device.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**2 - MASK COMPONENTS (section A)**

- a) cushion; b) frame; c) elbow connector; d) swivel fitting; e) headgear; f) buckle.

**3 - ASSEMBLING AND WEARING THE MASK (section B)**

- Bring the gray cushion (a) ring closer to the frame (b) so as to first couple the upper part of the cushion to the upper tooth of the frame, then place your fingers and push to also hook the lower part of the cushion to the frame (B.1).
- Holding the mask gently in close contact to the face, wear the headgear (e) over the head making sure the cushion (a) adheres to the profile around the nose (B.2).
- Attach the two lower straps of the headgear (e) to the hooks of the frame (b) so that they pass under the ears (B.3).
- Slightly pull the straps to ensure the stability of the mask according to the sequence B.4.
- Connect the swivel fitting (d) to the equipment tube (B.5), then turn on the equipment at the prescribed pressure.
- Lie down and, if necessary, adjust the headgear again (B.6).

**4 - USING THE MASK**

To properly remove, disassemble and reassemble the mask, refer to the respective illustrations in section C - D - E.

**5 - CLEANING THE MASK**

For daily cleaning, the mask must be disassembled in its components, each of which (frame, cushion, elbow connector and its exhalation orifices) must be carefully washed in warm water (about 30°C) using neutral soap and then rinsed thoroughly with clean fresh water. Finally, the parts must be air dried, away from direct sunlight.

If necessary, carry out these operations with the help of a soft-bristled toothbrush.

**To clean the headgear, perform the same procedure weekly.**

For a deeper decontamination of the mask, consult the "Guide to decontamination, disinfection and sterilization" available electronically on the website [it.medicaldevice.airliquide.com](http://it.medicaldevice.airliquide.com) or at your dealer, together with the instructions provided by the manufacturer of the disinfectant solution, especially regarding the compatibility with the materials of the various components.

**⚠ PRECAUTIONS**

- Remove the mask after disconnecting at least one of the lower headgear straps.
- Do not open the straps to avoid having to adjust the headgear again the next time you use it.
- Be very careful to remove any dirt from the exhalation holes.
- The headgear can only be washed in water.
- Do not clean the mask with solutions containing alcohol, aromatic compounds, humidifiers, antibacterial agents, hydrating agents or aromatic oils as they would damage the product and reduce its life.
- Do not leave the mask immersed for more than 10 minutes.
- Do not clean the mask in the dishwasher.
- Check carefully all the components to make sure there are no soap deposits or foreign bodies inside the mask.

**6 - REGENERATION OF THE MASK FOR REUSE ON DIFFERENT PATIENTS**

The mask can be regenerated for up to 20 cycles by one of the following processes:

- high level thermal disinfection;
- high level chemical disinfection.

For full detailed mask regeneration instructions, consult the “Guide to decontamination, disinfection and sterilization” available electronically on the website [it.medicaldevice.airliquide.com](http://it.medicaldevice.airliquide.com) or at your dealer.

**⚠ PRECAUTIONS**

- *The mask must always be regenerated before use on a new patient.*
- *The headgear must be replaced.*

**7 - TROUBLESHOOTING**

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
The airflow doesn't reach the mask.	<ul style="list-style-type: none"> <li>- Flow generator turned off or not working.</li> <li>- Breathing circuit ventilator not inserted correctly.</li> </ul>	<ul style="list-style-type: none"> <li>- Switch on or adjust the flow generator.</li> <li>- Connect the breathing circuit ventilator to the flow generator and the mask.</li> </ul>

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
Noisy and annoying air leaks.	<ul style="list-style-type: none"> <li>- The mask has not been assembled correctly.</li> <li>- The mask is not positioned correctly.</li> <li>- The size of the mask may not be correct.</li> </ul>	<ul style="list-style-type: none"> <li>- Disassemble the mask and then reassemble it as described in the user manual.</li> <li>- Reposition the mask on your face and adjust the headgear straps as described in the user manual.</li> <li>- Contact your doctor or therapy technician to find the right size.</li> </ul>
The exhalation orifices become dirty to the point of occlusion.	Incorrect or insufficient daily cleaning.	<p>Follow the indicated procedure:</p> <ul style="list-style-type: none"> <li>- Disassemble the mask in all its components;</li> <li>- Leave the elbow connector with the exhalation orifices immersed in a bath of hot water (at maximum domestic water temperature) and neutral soap for about 2 hours;</li> <li>- Rinse the component with fresh clean water and pay particular attention when removing any soap deposits;</li> <li>- Dry it with a dry cloth/ allow to air dry, away from direct sunlight;</li> <li>- Assemble the mask.</li> </ul> <p>If necessary, carry out these operations with the aid of a soft bristle brush and brush the exhalation orifices for 30 seconds under running water.</p>


**8 - DISPOSAL**

The mask components are made of the materials indicated below. They do not contain hazardous substances and can therefore be disposed of with ordinary household waste.

COMPONENT	MATERIAL
Cushion (a)	Silicone + PBT
Frame (b)	Polycarbonate
Elbow connector (c)	Polycarbonate + PBT + Polyester
Swivel fitting (d)	Copolyester
Headgear (e)	Nylon/UBL/Polyurethane
Buckle (f)	Nylon

**9 - TECHNICAL DATA**

**PRESSURE-FLOW CURVE**

*\*Due to machining differences, the outflow from the exhalation orifices may vary (± 10%).*

<b>Regulatory information</b>	CE marking in accordance with Regulation 2017/745 (EU) and subsequent implementation amendments. Class IIa Main product standard UNI EN ISO 17510
<b>Operating pressure range</b>	4 cmH <sub>2</sub> O / 32 cmH <sub>2</sub> O
<b>Operating temperature range</b>	+5°C / +40°C
<b>Dead space</b>	Size S: 64 mL Size M: 80 mL Size L: 96 mL Size XL: 115 mL Size M Asia: 74 mL
<b>Resistance</b>	at 50 l/min: 0.15 cmH <sub>2</sub> O at 100 l/min: 0.82 cmH <sub>2</sub> O
<b>A-Weighted Sound Pressure (UNI EN ISO 3744)</b>	12.6 ± 3.2 dB(A)
<b>A-Weighted Sound Power (UNI EN ISO 3744)</b>	20.6 ± 3.2 dB(A)
<b>Dimensions (without headgear)</b>	Size S: 150 x 80 x 80 mm Size M: 150 x 80 x 85 mm Size L: 150 x 80 x 90 mm Size XL: 150 x 80 x 95 mm Size M Asia: 150 x 80 x 80 mm
<b>Weight</b>	Size S: 75 g Size M: 80 g Size L: 87 g Size XL: 95 g Size M Asia: 78 g

**10 - WARRANTY CERTIFICATE**

Air Liquide Medical Systems S.r.l. guarantees its product against material or manufacturing defects for 180 days from the date of purchase, provided that the conditions of use in the instructions for use are observed. If the product is faulty under normal use, Air Liquide Medical Systems S.r.l. will, at its discretion, repair or replace the faulty product or its components.

## ENGLISH

The manufacturer shall be responsible for the safety, reliability and performance of the device if it is used in accordance with the operating instructions and its intended use.

Keep the original casing for the entire warranty period. The warranty in case of product defects can be applied by contacting your dealer.

For more information on the rights under the warranty, please contact your retailer.

### COPYRIGHT

All the information in this manual may not be used for purposes other than the original. This manual is the property of Air Liquide Medical Systems S.r.l. and it may not be reproduced, in whole or in part, without written permission from the company. All rights are reserved.

### TECHNICAL CHARACTERISTICS UPDATES

For a continuous improvement of performance, safety and reliability, all medical devices produced by Air Liquide Medical Systems S.r.l. are regularly subjected to review and change.

Instruction manuals are therefore modified to ensure their constant consistency with the characteristics of devices placed on the market.

If the instruction manual accompanying this device is lost, contact the manufacturer to get a copy of the version corresponding to the device by quoting the references shown on the label.

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AIR LIQUIDE MEDICAL SYSTEMS S.r.l.