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## ENGLISH

### INDICATION FOR USE & INTENDED USERS

Respiro Hospital F2 Non Vented Disposable is an oronasal mask without calibrated exhalation orifices intended to be used, for short periods of time (less than 7 days) in a hospital setting, by a single adult patient (weighing over 30kg) in controlled non-invasive positive pressure ventilation (NIPPV).

### WARNINGS

- The mask should only be used with a ventilator fitted with an active exhalation valve and/or with a respiratory circuit with an exhalation valve recommended by your doctor or respiratory therapist.
- In case of discomfort, irritation or allergic reactions to any component of the mask, consult your doctor or respiratory physiotherapist.
- When administering additional oxygen, smoking or open flames are prohibited.
- Turn off the oxygen dispenser when not using the equipment.
- Do not use the mask in case of vomiting or nausea.
- Keep away from light.
- Before use, check that the mask is in good condition. If damaged during transport, notify your dealer.
- Do not leave the components of the mask unattended, some of these may be swallowed by children.
- The device is a single use device; reuse on multiple patients may result in transmission of infections, serious harm or even death of the patient.

Contact your dealer in their capacity as a healthcare professional if there are issues in terms of operation or size, or if you have questions regarding the application of the medical device.

Should a serious accident occur in connection with the device, the manufacturer and the competent authority of the Member State in which the user and/or patient is established must be notified.

## 2. TROUBLESHOOTING

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
There is no airflow in the mask.	- Flow generator turned off or not working. - Breathing circuit ventilator not inserted correctly.	- Switch on or adjust the flow generator.  - Connect the breathing circuit ventilator to the flow generator and the mask.
There are noisy and annoying air leaks.	- The mask has not been assembled correctly.  - The mask is not positioned correctly.  - The size of the mask may not be correct.	- Disassemble the mask and then reassemble it as described in the user manual. - Reposition the mask on your face and adjust the headgear straps as described in the user manual. - Contact your doctor or therapy technician to find the right size.

## 3. DISPOSAL

The mask components are made of the materials indicated below. They do not contain hazardous substances and can therefore be disposed of with ordinary household waste.

COMPONENT	MATERIAL
Cushion (a)	Silicone
Frame (b)	Polycarbonate
Forehead support (c)	Silicone
Blue elbow connector (d)	Polycarbonate
Headgear (e)	Nylon/UBL/Polyurethane
Clip (f)	Nylon

## 4. TECHNICAL DATA

Regulatory information	CE marking in accordance with Regulation 2017/745 (EU) as amended and updated. Class IIa Main product standard UNI EN ISO 5356-1
Operating pressure range	4 cmH <sub>2</sub> O / 32 cmH <sub>2</sub> O
Operating temperature range	+5 °C / +40 °C
Dead space	S: 284 ml - M: 294 ml - L: 325 ml
Dimensions (without headgear)	S: 185x100x110 mm - M: 185x105x110 mm - L: 195x100x110 mm
Weight	S: 129 g - M: 132 g - L: 133 g

## 5. WARRANTY CERTIFICATE

Air Liquide Medical Systems S.r.l. guarantees its product against material or manufacturing defects for 180 days from the date of purchase, provided that the conditions of use in the instructions for use are observed. If the product is faulty under normal use, Air Liquide Medical Systems S.r.l. will, at its discretion, repair or replace the faulty product or its components. The manufacturer shall be responsible for the safety, reliability and performance of the device if it is used in accordance with the operating instructions and its intended use. Keep the original casing for the entire warranty period. In the event of product defects, the warranty can be enforced by contacting your dealer. For more information on our rights under the warranty, please contact your dealer.

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