
MBS Basic

MBS Adult Set

MBS Junior Set

Air Liquide
HEALTHCARE

Via Bisceglie, 66 - 20152 Milano - Italy
Plant: Via Dei Prati, 62 - 25073 BOVEZZO (BS) (Italy)
Tel. 030.20159.11 - Fax 030.200051
medic@airliquide.com
airliquide.com

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Please read the instructions for use and the important notes (section 1.4) carefully before using the device.

1 - IMPORTANT INDICATIONS

1.1 - INTENDED USE

This medical device belongs to the category "accessories for aerosol therapy" and is intended to be used for administering a solution in aerosol form.

1.2 - INDICATIONS FOR USE

The device is a nebuliser intended to be used for the administration of drugs in aerosol form. The nebuliser must be connected to a pneumatic aerosol therapy device that generates compressed air. It can be used to treat diseases of the lower

airways (including bronchial asthma, chronic bronchitis, pulmonary emphysema and cystic fibrosis).

The administration of drugs requires a medical prescription specifying the type of drug, the dosage and the length of treatment.

MB5 is designed and made in Italy in accordance with EU Regulation 2017/745 on Medical Devices (as amended and updated) and bears the CE marking.

1.3 - SYMBOLS USED

	Important, read the instructions carefully		Attention!
	Lot number		Product code
	UDI code		Medical device
	Storage temperature		Production date
	Storage humidity		Single-patient use device
	Storage pressure		Manufacturer
	Does not contain phthalates		No traces of natural latex

1.4 - IMPORTANT NOTES

- Read the instructions for use carefully and keep them in a safe place for future reference.
- The nebuliser must be connected to a pneumatic aerosol therapy device that generates compressed air. Use only CE-marked devices.
- Single-patient use nebuliser, non-heat sterilisable. Must be used for one person only.
- The device must only be used for treating the lower airways in accordance with these instructions for use. It must not be used for any other purpose.
- There are no contraindications to the use of the device.
- Some components are small enough to swallow and may lead to suffocation.
- Read the list of components in the table in section 7 to see if there are any components that have caused an allergic reaction in the past.
- Consult the manufacturer's instructions of the device with which the nebuliser is used, to make sure they are compatible (see 6 - TECHNICAL DATA).
- Children must only be allowed to use the device under the supervision of an adult who has read and understands these instructions.
- The administration of drugs in aerosol form requires a medical prescription, specifying the type of drug, the dosage and the length of treatment.
- Do not use the nebuliser in the bath or shower.
- The manufacturer has defined a service life of at least 3 years from first use.
- Should a serious accident occur in connection with the device, the manufacturer and the competent authority of the Member State in which the user and/or patient is established must be notified.
- In case of doubt about the use of the product, contact the manufacturer or physician or pharmacist.

1.5 - PATIENT GROUPS AND ENVIRONMENTS OF USE

MB5 provides highly effective inhalation therapy for patients of all ages, from young children (used with paediatric mask) to cooperating adults (used with mouthpiece, nosepiece, or adult mask), and can be used in both home and hospital settings.

1.6 - DESCRIPTION OF THE DEVICE

MB5 is a nebuliser specifically for treatment of the lower airways. The nebuliser is for single-patient use, as it is in material that cannot be heat sterilised. The nebuliser must be connected to a pneumatic aerosol therapy device that generates compressed air.

The nebuliser is available in multiple configurations:

- MB5 Basic: nebuliser with mouthpiece
- MB5 Junior set: nebuliser with mouthpiece, paediatric mask with elastic band, connecting tube and connector
- MB5 Adult set: nebuliser with mouthpiece, adult mask with elastic band, connecting tube and connector

1.7 - CONTENTS

The package contains (Fig. A) the parts listed below. A nebuliser comprising the following:

- upper section (1)
 - atomiser (2)
 - lower section (3)
 - one mouthpiece (4)
 - instructions for use
 - one mask (5), if provided by the configuration. The mask is fitted with an elastic band for securing to the face.
 - one connecting tube with connector, if provided by the configuration.
- The lower section (3) is fitted with a terminal to which the air supply tube is connected.

2 - PREPARING AND USING THE MB5

MB5 requires connection to a source of compressed air. In all cases, it is important to make sure that the source of air is compatible with the nebuliser (see section 6 - TECHNICAL DATA).

Please read section 3 before using the device.

Proceed as follows to use the MB5:

- make sure the atomiser (2) is correctly inserted in the central sleeve in the lower section of the nebuliser (3).
- the lower section of the nebuliser (3) has a graduated scale that shows the approximate amount of liquid inside. The purpose is merely to give an approximate indication of the contents. Use a graduated syringe to get an accurate dose. Pour liquid into the tank (Fig. B) according to the amount recommended by your physician

Readings must be taken on the graduated scale with the atomiser mounted on the lower section of the nebuliser.

- tighten the upper section (1) of the nebuliser onto the lower section (3) (Fig.C)
- Insert the mouthpiece (4) to the upper section of the nebuliser (1) (Fig. A)
- if using the mask, the latter must be secured to the upper section of the nebuliser (Fig.D)
- connect the tube to the lower end of the nebuliser, and connect the other end to the air intake of the aerosol device. If necessary, use the black connector at the end of the tube so as to connect it to the aerosol device.
- start up the device.
- during treatment, hold the nebuliser upright in your hand to enable it to work efficiently
- if the above operations have been performed correctly, nebulised drug will come out of the nebuliser.

The mouthpiece (4) must be placed between closed lips, even during exhalation: exhaled air is released from a special hole at its base.

If using with mask (5), the latter must be gently placed on the face. There are openings on the side that allow exhaled air to be released. When you have used up all the drug or interrupt treatment, switch off the aerosol device or turn off the air supply.

A small amount of liquid may remain in the nebuliser at the end of the session but this has no importance.

After use, disconnect the nebuliser from the tube and the device from the mains.

3 - CLEANING AND DISINFECTING

Before proceeding, unplug the device from the power supply and disconnect the MB5 and tube from the device.

If you need to clean or disinfect this medical device using a method other than that given in this section, contact the manufacturer first for a copy of the data sheet.

To avoid the risk of microbial contamination, the MB5 must be cleaned and disinfected after each treatment and before use, following the instructions carefully. Before cleaning or disinfecting, components must be disassembled as indicated below:

- remove the tube
- remove the mouthpiece (4) or mask (5) from the nebuliser (1, 2, 3)
- undo the upper section of the nebuliser (1)
- remove the atomiser (2)

3.1 - CLEANING

The disassembled nebuliser (1, 2, 3), mouthpiece and mask (4, 4.a, 5) may be washed in hot water at an approximate temperature of 40°C with the addition of detergent (refer to the detergent manufacturer's instructions), then rinsed in running water.

When all the parts have been cleaned, dry them with a soft cloth and place in the storage compartment.

3.2 - DISINFECTING

After washing, all components of the nebuliser (1, 2, 3), the mouthpiece and mask (4, 5) may be immersed in a cold disinfectant solution (e.g. a solution containing up to 2% sodium hypochlorite). Follow the manufacturer's instructions.

4 - MAINTENANCE

When cleaning and disinfecting MB5 components, you may check them visually and determine when they need to be replaced. We recommend inspecting the appearance and condition of the components after each cleaning and disinfecting cycle.

The MB5 must be replaced when there are clear signs of deterioration.

5 - PROBLEMS, POSSIBLE CAUSES AND SOLUTIONS

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
The nebuliser does not generate a spray.	- The device is switched off.	- Switch on the device.
	- The liquid has not been poured into the nebuliser.	- Pour liquid into the nebuliser.
	- Too much liquid has been poured into the nebuliser.	- Reduce the quantity of drug in the nebuliser to below the maximum level.
	- The nebuliser has not been assembled correctly.	- Check that the atomiser (2) is properly inserted in the base of the nebuliser (3). Screw the upper section (1) onto the base and check that air comes out of the little hole at the top of the base.
	- The air outlet in the base (3) is blocked.	- Unblock the hole.

	- The air tube is not connected correctly	- Connect the air tube properly
The nebuliser does not perform as stated.	- The air source does not provide the required pressure and rate of flow	- Check the characteristics of the air source
The nebulised liquid does not reach the mask or there is difficulty exhaling	- The mouthpiece or nosepiece have not been fitted correctly	- Assemble the mask or mouth-piece correctly

6 - TECHNICAL DATA

Minimum supply pressure	1,5 bar	
Maximum supply pressure	3,2 bar	
Air flow at minimum pressure	7 l/min	
Air flow at maximum pressure	12 l/min	
Maximum filling volume	8 ml	
Recommended filling volume	3 ml	
MMAD, Diagram 1 (*) (**)	4,45 µm	
GSD (*) (**)	2,80 µm	
Breathable fraction (%<5 µm) (*) (**)	54,88%	
Particle distribution (*) (**)		
%>5 µm	45,12%	
2 µm <%< 5 µm	29,54%	
%< 2 µm	25,34%	
Residual volume	1,6 ml	
	At minimum pressure and flow	At maximum pressure and flow
Aerosol Output	1,3 ml	1,4 ml
Aerosol Output Rate	0,16 ml/min	0,30 ml/min
Percentage of filling volume per minute	5,2%	10%

(*) The following information on nebuliser performance is based on tests using adult ventilator schemes and is likely to differ from those claimed for paediatric or infant populations. The use of a high-viscosity suspension and/or solution may alter the particle size distribution curve, the MMAD, aerosol emission and/or aerosol emission rate, which may therefore differ from those disclosed by the manufacturer. Tests performed with 0,9% NaCl saline solution (**). The MMAD values shown are measured with the Malvern Spraytec instrument and using NaCl 0,9% saline.

7 - DISPOSAL OF THE MB5

If taken out of use, the device and/or its components must be disposed of in accordance with the regulations in force and environmental protection. The table shows the different types of device components. Check the regulations in force in your municipality. None of the materials used contains phthalates or traces of natural Latex.

COMPONENT	MATERIAL
Upper section (1)	Polypropylene
Atomiser (2)	Polypropylene
Lower section (3)	Polystyrene
Mouthpiece (4)	Polypropylene
Mask (5)	Non-toxic PVC
Tube	Non-toxic PVC
Connector	Nylon

8 - OPERATING CONDITIONS

- temperature from +5°C to +40°C
- relative humidity from 15% to 93%
- atmospheric pressure from 700 hPa to 1060 hPa

9 - STORAGE CONDITIONS

- temperature from -25°C to +50°C
- relative humidity less than 93%
- atmospheric pressure from 700 hPa to 1060 hPa

10 - WARRANTY CONDITIONS

The manufacturer is liable for the safety, reliability and performance of the device, provided that it is used in accordance with the instructions and for the intended use only, and any repairs are carried out by the manufacturer or an approved service centre. The manufacturer declines all liability and the warranty will be invalidated if the user fails to follow the instructions and the important notes contained herein.

11 - TECHNICAL UPDATES

Air Liquide Medical Systems S.r.l. periodically reviews and modifies all its medical devices in order to improve their performance, safety and reliability. The instruction booklets are updated accordingly to include any new or changed features. If the booklet accompanying the device gets damaged or mislaid, a replacement copy can be obtained from the manufacturer by quoting the data shown on the rating plate.

12 - COPYRIGHT

The information contained in this booklet must not be used for any purpose other than the reason for which it is provided. This booklet is the property of Air Liquide Medical Systems S.r.l. and no part of it may be reproduced without the manufacturer's written permission. All rights reserved.