NEBULA smart

ISTRUZIONI PER L'USO ● INSTRUCTIONS FOR USE ● MODE D'EMPLOI ● BEDIENUNGSANLEITUNG ● BRUKSANVISNING

Air Liquide HEALTHCARE

ITALIANO ● ENGLISH ● FRANÇAIS ● DEUTSCH ● NORSK




Fig. A

11
6
13
12
20
14
22
9
19
10
18
15
16
17
NON SMONTARE QUESTE DUE PARTI



Fig. B


Fig. C1

Fig. C2


Fig. D


Fig. E


Fig. F



Fig. G

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⚠ Please read the instructions for use and the important notes (paragraph 1.3) carefully before using the device.

1 - IMPORTANT INDICATIONS

1.1 - INTENDED USE

This medical device belongs to the category “aerosol therapy devices” and is intended to be used for administering a solution in aerosol form.

1.2 - INDICATIONS FOR USE AND INTENDED USERS

The device is an aerosol therapy device intended to be used for the administration of drugs in aerosol form.
Used with the nebuliser, the device can be used to treat diseases of the lower airways (including bronchial asthma, chronic bronchitis, pulmonary emphysema and cystic fibrosis).
The micronized nasal douche is used for treating infections of the upper airways and for physiological flushing (with saline, hypertonic solution, or thermal spring water) and/or therapeutic flushing (with drugs) of the nasal cavities.

⚠ The administration of drugs requires a medical prescription specifying the type of drug, the dosage and the length of treatment.

This aerosol therapy device and its accessories have been designed and made in Italy in accordance with Regulation (EU) 2017/745 regarding medical devices (and subsequent updates) and bear the CE marking.

1.3 - IMPORTANT NOTES

- Read the instructions for use carefully and keep them in a safe place for future reference.
- Do not use the device and/or accessories for anything other than the intended use.
- There are no contraindications to the use of the device.
- The administration of drugs requires a medical prescription specifying the type of drug, the dosage and the length of treatment.
- Children must only be allowed to use the device under the supervision of an adult who has read and understands these instructions.
- Some components are small enough to swallow and may lead to suffocation.
- The power cable and connecting tube are long enough to pose a strangulation hazard.
- Read the list of components in the table (section 7) to see if there are any components that have caused an allergic reaction in the past.
- Do not use essential oil solutions (e.g. menthol, eucalyptus, etc.) with Rinowash as they are not compatible with the material of which the device is made.
- During operation, the device must stand on a firm flat surface.
- During operation, never cover the air vents.
- During operation, keep the device well away from sources of heat.
- Do not touch the plug or the electrical components of the device with wet or damp hands.
- Device protected against finger access to hazardous parts and protected against vertically falling drops (IP21).
- Do not pour water or other liquids onto the device.
- Do not use the device in the bath or shower.
- If the device should fall into water, disconnect the power plug before removing it. After such an event, the device may not be used and requires a complete overhaul.
- Beware of compressed air coming out of the air outlet as this may cause a hazard; it must only be used to supply the nebuliser.
- The device is not suitable for use in the presence of anaesthetic mixtures with air, oxygen or nitrous oxide.
- The device is not suitable for use with a ventilator.
- Performance of the device is only guaranteed if all the accessories used are original ones. Only original accessories must be used.
- The manufacturer has defined a service life of at least 5 years from first use.
- Should a serious accident occur in connection with the device, the manufacturer and the competent authority of the Member State in which the user and/or patient is established must be notified.
- In case of doubt about the use of the product, contact the manufacturer or physician or pharmacist.

1.4 - PATIENT GROUPS AND ENVIRONMENTS OF USE

The NEBULA SMART aerosol therapy device provides highly effective inhalation therapy for patients of all ages, from small children to adults, and can be used in both home and hospital settings.

1.5 - DESCRIPTION OF THE DEVICE

The NEBULA SMART aerosol therapy device is made up of a compressor driven by an electric motor in the basic unit, a SoCu nebuliser (with mouthpiece, PVC mask for adults and children, nosepiece), a Rinowash micronized nasal douche, and a flexible tube connecting the nebuliser/nasal douche to the basic unit.

1.6 - CONTENTS

The package contains (Fig. A): The basic unit (1), ON/OFF switch (2), air intake (3), nebuliser support (4), power cable plug (5), filter slot (6).
It also contains:

- SoCu nebuliser made up of the following:
 - nebuliser (upper section (7), lower section (8))
 - mouthpiece with exhalation valve (14)
 - mask for adults (10), paediatric mask (9)
 - nosepiece (22)
- connecting tube (11)
- spare air filters (12)
- power supply cable (13)
- one Rinowash micronized nasal douche made up of the following:
 - one base with button (15), connecting element for the connecting tube (16), base cap (17), and central tank (18)
 - atomiser (19)
 - outer bell (20)
- Protective cover (21)

2 - PREPARATION AND USE

2.1 - PREPARING THE BASIC UNIT

⚠ Before use, make sure the voltage indicated on the rating plate underneath the device corresponds to the mains voltage, and the wiring system complies with current regulations.
Do not use the device if the cable is damaged.
Check that the casing is intact and that there is no sign of dents or damage.

Take the power cable (13) and plug it into the power outlet (5).
Make sure the switch (2) is in position “O” (OFF). Connect the cable to the main unit and the electric plug in the mains socket, press position “I” (ON) on the switch, and check that air comes out of the air intake (3). After checking, press position “O” (OFF) on the switch.

2.2 - USING THE SOCU

Please read paragraph 3.2 before using the device.
Proceed as follows to use the SoCu nebuliser:

- The lower section of the nebuliser (8) has a graduated scale that shows the approximate amount of liquid inside. The purpose is merely to give an approximate indication of the contents. Use a graduated syringe to get an accurate dose. Pour liquid into the tank (Fig. B) according to the amount recommended by your physician
- Attach the upper section of the nebuliser (7) to the base (8) and turn it until it clicks (Fig. C1).
- Connect the lower end of the tube (11) to the nebuliser, and the other end to the air intake (3) of the device
- Connect the mask (9/10), mouthpiece (14), or nosepiece (22) to the base of the nebuliser (Fig. C2)
- Press position “I” (ON) on the switch (2)
- During treatment, hold the nebuliser upright in your hand to enable it to work efficiently
- If the above operations have been performed correctly, nebulised liquid will come out of the mask

When you have used up all the liquid or wish to temporarily interrupt treatment, press “O” (OFF) on the switch and place the nebuliser back on its support (4) (Fig. A).
A small amount of liquid may remain in the nebuliser at the end of the session but this has no importance.
If a small amount of moisture remains inside the air hose at the end of the session, leave the device switched on and allow it to dry by ventilation.
The presence of liquid inside the tube is the result of water vapour forming under certain environmental conditions.

2.3 - USING THE RINOWASH MICRONIZED NASAL DOUCHE

Please read paragraph 3.2 before using the device.
Proceed as follows to use the Rinowash micronized nasal douche:

- Insert the liquid (maximum 15 ml) into the tank (Fig. D) checking that the liquid matches the graduated scale.

⚠ Readings must only be taken on the graduated scale when the atomiser is not mounted into the tank

- Fit the atomiser (19) onto the sleeve of the tank. Make sure it clicks into position
- Fit on the outer bell (20), turning slightly, to cover the seal
- Attach one end of the connecting tube (11) to the Rinowash connector (16) and the other end to the air intake (3)
- Press “I” (ON) on the switch (2) to start up the device
- Bring the Rinowash to the nose using the hand opposite the nostril to be treated, and rest the end against one nostril (Fig. E). Keep your head tilted forward.
- Press the button on the base (15) to generate a spray
- Therapeutic spray is delivered through the hole in the outer bell (20), and residual fluid collects in the chamber between the outer bell and the tank (Fig. F)
- During treatment, breathe normally through the nose. If mucus forms in the nostril, blow it out. The space between the outer bell and the tank acts as a chamber for collecting the fluid that drains out of the nose.
- Treat the nostrils alternately, as required

Nebulisation is very rapid and the session only lasts a few minutes.

⚠ Make sure all the parts are connected firmly together. Check that the ends of the tubes are inserted properly and the components of the Rinowash are mounted correctly.
If the device is not assembled correctly, micronization of the liquid may be affected. If this happens, rotate the atomiser (19).

As stated under Important Notes, some liquids are not compatible with the materials of which the nebuliser, Rinowash, and spacer mouth mask are made. They may interact and cause damage (the materials are shown in section “7 - DISPOSING OF THE DEVICE AND/OR ITS COMPONENTS”).

3 - CLEANING, DISINFECTING AND STERILISING

⚠ Before proceeding, unplug the device and disconnect the nebuliser, Rinowash nasal douche, and tube

⚠ If it is necessary to clean, disinfect, or sterilise this medical device using methods other than those indicated below in Chapter 3, verify that this is possible by consulting the technical data sheet available from the manufacturer.

3.1 - CLEANING, DISINFECTING AND STERILISING THE BASIC UNIT

3.1.1 - CLEANING THE BASIC UNIT

The basic unit (1) must not be cleaned in water as it could get inside and damage the electrical components, which would cause a safety hazard.
Clean the basic unit (1) with a soft damp cloth, with a small amount of mild (non-abrasive) detergent, if necessary.

3.1.2 - DISINFECTING THE BASIC UNIT

The basic unit can be disinfected using a soft sponge or cloth and some cold disinfectant solution (e.g. a solution containing up to 2% sodium hypochlorite). Follow the manufacturer’s instructions.
Make sure the solution is compatible with the material of which the basic unit is made (warning: consult the important notes verifying that the disinfectant solution is compatible with the materials used, given in section “7 - DISPOSING OF THE DEVICE AND/OR ITS COMPONENTS”).

3.1.3 - STERILISING THE BASIC UNIT

The basic unit (1) must not be sterilised.

3.2 - CLEANING, DISINFECTING AND STERILISING THE ACCESSORIES

To avoid the risk of microbial contamination, accessories must be cleaned, disinfected, and, if necessary, sterilised after each treatment and before use, carefully following the instructions.
Before cleaning, disinfecting or sterilising any of the components, disassemble them as follows:

- remove the tube (11)
- remove the accessories from the nebuliser
- undo the upper section of the nebuliser (7)

For the Rinowash:

- remove the tube (11)
- remove the outer bell (20) by turning gently
- remove the atomiser (19) by pressing on the tab (Fig. G)
- remove the base cap (17).
- The seal at the base of the tank (18) must not be removed.
- The main section and the tank (18) form a single component and cannot be disassembled.

3.2.1 - CLEANING THE ACCESSORIES

The disassembled accessories can be washed in water at an approximate temperature of 40°C with the addition of detergent (refer to the manufacturer's instructions), then rinsed in running water.
When all parts have been cleaned, they can be dried with a soft cloth and stored.

3.2.2 - DISINFECTING THE ACCESSORIES

After being washed, the accessories can be immersed in a cold disinfectant solution (e.g. a solution containing up to 2% sodium hypochlorite) following the instructions given by the manufacturer.
The nebuliser and its accessories (except the air hose) can be disinfected with chemical disinfectants (e.g. a solution containing up to 2% sodium hypochlorite following dosage and use limitations as provided by the disinfectant manufacturer).
The nebuliser and its accessories (except hose and masks) can be disinfected by boiling them in water for a maximum time of 10 min: make sure that the components are not in direct contact with the bottom of the pot.

3.2.3 - STERILISING THE ACCESSORIES

The nebuliser and its accessories (except air hose and masks) can be heat steam sterilised up to 121°C for 18 minutes for a maximum of 20 cycles. After sterilisation, always let all components cool down to ambient temperature before further use. Do not repeat the sterilisation cycle while components are still warm.
The sterilisation methods have been validated and documented by the manufacturer in Technical Test S-4/12 no. 317. Repeated sterilisation cycles do not affect the quality and performance of the product. The manufacturer does not guarantee the sterile state of the device.

3.3 - CLEANING, DISINFECTING AND STERILISING THE TUBE

The tube (11) needs to be cleaned after each session. Detach it from the nebuliser and leave the device on until all signs of moisture have been eliminated. If the tube is particularly dirty, it needs to be replaced.

⚠ The tube (11) must not be sterilised hot or in an autoclave or boiled.

4 - MAINTENANCE

The NEBULA aerosol therapy device does not require any particular maintenance. The compressor does not need to be lubricated as the piston is the dry-operating type.
A special external filter keeps out dust and moisture, and should be replaced periodically.

⚠ Always unplug from the mains before carrying out the following operations.

4.1 - REPLACING THE AIR FILTER

The air filter needs to be replaced if it gets dirty, and once a year if the device is used frequently.
To replace the air filter (6), remove the cap using a screwdriver and insert a new one. The filter cannot be washed or regenerated. It must always be replaced with a new one.

4.2 - MAINTENANCE OF THE ACCESSORIES

When cleaning, disinfecting or sterilising the components, you can check them visually and determine when they need to be replaced. Repeated sterilisation may affect the materials of which the accessories are made; after each sterilisation cycle, inspect the appearance and integrity of the components.
The accessories must be replaced when there are clear signs of deterioration. In any case, we recommend that accessories be replaced at least once a year or, in the case of frequent use, replaced as their features change.

GLISH • ENGLISH • ENGLISH • ENGLISH • ENGLISH • ENGLISH

5 - PROBLEMS, POSSIBLE CAUSES AND SOLUTIONS

TROUBLE-SHOOTING CHART

PROBLEMS	POSSIBLE CAUSES	POSSIBLE SOLUTIONS
The device does not work	- The plug is not inserted in the socket	- Insert the plug in the socket and press position “I” (ON) on the switch (2)
	- The switch (2) is in position “O” (OFF)	- Press position “I” (ON) on the switch (2)
	- The thermostat has tripped	- Check that the air vents are not clogged
The nebuliser does not generate a spray	- Device off	- Switch on the device
	- The liquid has not been poured into the nebuliser	- Pour liquid into the nebuliser
	- Too much liquid has been poured into the nebuliser	- Reduce the quantity of liquid in the nebuliser to below the maximum level
	- The nebuliser has not been assembled correctly	- Screw the upper section (7) onto the base and check that air comes out of the little hole at the top of the nebuliser base
	- The small hole in the air outlet in the base (8) is clogged	- Unclog the small hole
	- The tube (11) is not connected securely	- Connect the tube (11) securely to the air intake on the device (3) and the air intake on the nebuliser
The Rinowash does not deliver a spray	- Device off	- Switch on the device
	- There is no liquid in the tank	- Pour liquid into the tank
	- The tube is not connected securely or the Rinowash has not been assembled correctly	- Connect the tube securely to the air intake on the device and the air intake on the Rinowash. Slightly turn the atomiser. Make sure all the parts have been assembled correctly.
Liquid leaks from the base of the bell	- The outer bell has not been fitted on properly and the seal does not hold	- Fit on the outer bell (20) correctly, turning slightly. Make sure the whole seal is covered


If you have failed to solve the problem after performing the above checks and the device still does not work, have it checked by an authorised service centre.

6 - TECHNICAL DATA


Supply voltage	230 V ~ 50 Hz
Operation	Continuous
Nominal air pressure to the device	2.5 bar
Nominal air flow to the device	10 l/min
Maximum input	126.5 VA
Safety fuse	F 1.6 A - 250 V
Dimensions of the basic unit	230x180x185H mm
Weight of the basic unit	3 kg
Noise level at 1 metre	60 dB
Recommended nebuliser fill volume	3 ml
Maximum nebuliser fill volume	9 ml

Minimum Rinowash fill volume	1.5 ml
Maximum Rinowash fill volume	15 ml
Aerosol output with nebuliser (*)	0.49 mg
Aerosol output rate with nebuliser (*)	0.11 mg/min
MMAD with nebuliser, graph 1 (*)	3.95 µm
GSD with nebuliser (*)	2.20 µm
Breathable fraction (%<5 µm) (*)	60.2%
Particle distribution (*)	
%>5 µm	37.5%
2 µm <%< 5 µm	41.1%
%< 2 µm	19.2%
Percentage of fill volume emitted per minute (*)	3.6%
Residual volume:	1.38 mg
MMAD with Rinowash micronized nasal douche, graph 2 (*)(**)	49 µm (**)

GRAPH 1


GRAPH 2

(*) The following information on nebuliser performance is based on tests using adult ventilator schemes and is likely to differ from those claimed for paediatric or infant populations. The use of a high-viscosity suspension and/or solution may alter the particle size distribution curve, the MMAD, aerosol emission and/or aerosol emission rate, which may therefore differ from those disclosed by the manufacturer.

(**) The MMAD values shown are measured with the Malvern Spraytec instrument and using NaCl 0.9% saline.

6.1 - SYMBOLS USED

	Class II device		Alternating current
	Device with type BF applied part	I	Device on
SN	Serial Number	O	Device off
UDI	UDI code		Manufacturer
REF	Product code		Production date
MD	Medical device		No traces of natural latex
	Warning!		Does not contain phthalates
	Warning: read the instructions carefully (blue symbol)		Storage temperature
	Keep the device dry		Storage humidity
IP21	Device protected against finger access to hazardous parts and protected against vertically falling drops		Storage pressure
	The device must be dismantled at the end of its lifecycle and disposed of in special bins. Do not dispose of in ordinary town waste bins (European Directive 2012/19/EU)		

7 - DISPOSING OF THE DEVICE AND/OR ITS COMPONENTS

If taken out of use, the device and/or its components must be disposed of in accordance with the regulations in force and environmental protection. The table shows the different types of device components.
Check the regulations in force in your municipality. Do not dispose of WEEE waste with regular household waste.
None of the materials used contains phthalates or traces of natural Latex.

COMPONENT	MATERIAL
Basic unit	Waste from Electric and Electronic Equipment (WEEE)
External casing of the basic unit	ABS
Cover	ABS
Nebuliser	Polypropylene
Adult/paediatric mask	Non-toxic PVC
Mouthpiece / Nosepiece	Polypropylene
Exhalation valve	Silicone
Rinowash body	Polycarbonate (some metal parts)
Rinowash base cap	Elastomer
Rinowash atomiser	Polycarbonate
Rinowash outer bell	Polycarbonate and silicone
Connecting tube	Non-toxic PVC (polyvinyl chloride)
Power cable	Waste from Electric and Electronic Equipment (WEEE)
Box	Cardboard
Instructions	Paper
Protective bag	Low-density polyethylene (LDPE)

8 - ELECTROMAGNETIC COMPATIBILITY

This device complies with the requirements of CEI EN 60601-1-2:2015 on the electromagnetic compatibility of medical devices. Compliance with standards on electromagnetic compatibility does not ensure complete product immunity; operation of this device may be affected by electromagnetic interference above the levels specified in IEC 60601-1-2. Contact the manufacturer for more information.

9 - OPERATING AND STORAGE CONDITIONS

9.1 - OPERATING CONDITIONS

- temperature from +5°C to +40°C
- relative humidity from 15% to 93%
- atmospheric pressure from 700 hPa to 1060 hPa

9.2 - STORAGE CONDITIONS

- temperature from -25°C to +50°C
- relative humidity less than 93%
- atmospheric pressure from 700 hPa to 1060 hPa

10 - SPARE PARTS

Use only original spare parts.

11 - WARRANTY CONDITIONS

The manufacturer is liable for the safety, reliability and performance of the device, provided that it is used in accordance with the instructions and for the intended use only, and any repairs are carried out by the manufacturer or an approved service centre. The manufacturer declines all liability and the warranty will be invalidated if the user fails to follow the instructions and the important notes contained herein.

12 - TECHNICAL FEATURES UPDATE

Air Liquide Medical Systems S.r.l. periodically reviews and modifies all its medical devices in order to improve their performance, safety and reliability. The instruction booklets are updated accordingly to include any new or changed features. If the booklet accompanying the device gets damaged or mislaid, a replacement copy can be obtained from the manufacturer by quoting the data shown on the rating plate.

13 - COPYRIGHT

The information contained in this booklet must not be used for any purpose other than the reason for which it is provided.
This booklet is the property of Air Liquide Medical Systems S.r.l. and no part of it may be reproduced without the manufacturer's written permission. All rights reserved.

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